

FILAPROP IFU

(62 X 300) mm

Date : 06.06.2024

<div>INSTRUCTIONS FOR USE</div> <div>FILAPROP™</div>	
<div>MANUFACTURED BY</div> <div><div>Meril</div><div>Endo-Surgery</div></div> <div><div><div><div></div><div>MERIL ENDO SURGERY PVT. LTD.</div></div><div>THIRD FLOOR, E1-E3 MERIL PARK, SURVEY No. 135/2/B &amp; 174/2 MUKTANAND MARG, CHALA, VAPI - 396 191, GUJARAT, INDIA</div><div>Customer Care No. : 18004194433</div><div>E-mail : enquiry.endosurgery@merilife.com</div><div>Web. : www.merilife.com</div></div><div><div><div>CE</div><div>0297</div></div><div><div>EC REP</div><div>OBELIS S.A</div></div><div>Bd. Général Wahis, 53, 1030 Brussels, Belgium</div><div>T : +32 2 732 59 54 F. : +32 2 732 60 03</div><div>E : mail@obelis.net W. : www.obelis.net</div></div><div>Document No. : ME/IFU/194</div><div>Rev.00 Rev. Dt. 2024/06</div></div>	
<div><div>FILAPROP™</div><div>NON-ABSORBABLE POLYPROPYLENE SURGICAL SUTURE U.S.P.</div></div> <div><div>DESCRIPTION</div><p>FILAPROP™ suture is a monofilament, synthetic non-absorbable, sterile surgical suture composed of an isotactic crystalline stereoisomer of polypropylene, a synthetic linear polyolefin. The molecular formula of polypropylene is (C<sub>3</sub>H<sub>5</sub>)<sub>n</sub>. FILAPROP™ sutures are either undyed or dyed with phthalocyanine blue.</p><p>FILAPROP™ is available in a range of gauge sizes and lengths attached to standard stainless steel needles of various types and sizes.</p><p>FILAPROP™ complies with the requirements of the United States Pharmacopeia for "Non Absorbable Surgical Suture" and the European Pharmacopeia for "Sterile Synthetic Non Absorbable monofilament Strand." However, it may be slightly oversize in diameter to U.S.P. for some suture sizes. U.S.P except for diameter.</p></div> <div><div>INDICATION</div><p>FILAPROP™ Sutures are indicated for approximation and/or ligation of:</p><ul style="list-style-type: none"><li>- Skin tissue</li><li>- Soft tissues including cardiac</li></ul></div> <div><div>INTENDED USER GROUP</div><p>FILAPROP™ Sutures are intended to be used by Healthcare professionals, general physicians and surgeons who performs skin closures, soft tissue approximation and/or ligation.</p></div> <div><div>INTENDED TARGET GROUP</div><p>FILAPROP™ Suture is intended to be used in both men and women and in all general populations including infants, children, and adults of all age groups.</p></div> <div><div>INTENDED PURPOSE</div><p>FILAPROP™ Sutures are intended to be used for skin closure and soft tissue (internal tissue) approximation and/or ligation.</p></div> <div><div>SELECTION CRITERIA</div><p>The selection of suture for implantation depends on patient condition, surgical experience, surgical technique and wound size.</p><p>The device is used by professional healthcare person only.</p></div> <div><div>PERFORMANCE</div><p>FILAPROP™ suture initially elicits a minimal inflammatory reaction in</p></div>	<div><p>leakage or resurgery.</p><p>4. Care should be taken to avoid damage while handling surgical needles. Grasp the needle in an area one third (1/3) to one half (1/2) of the distance from the attachment end to the sharp point. Grasping in the sharp point area could impair the penetration performance and cause fracture of the needle. Grasping at the butt or attachment end could cause needle-suture detachment, bend or breakage.</p><p>5. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.</p><p>6. Users should exercise caution while handling surgical needles to avoid inadvertent needle stick injury.</p><p>7. Discard the used needles appropriately.</p><p>8. Broken needles may result in foreign body residual which may require additional surgery.</p><p>9. Reuse of needle may result in transmission of blood borne pathogens to patients and users.</p></div> <div><div>ADVERSE REACTIONS</div><p>Adverse reactions associated with the use of FILAPROP™ include</p><ul style="list-style-type: none"><li>• Transitory local irritation at the wound site or transitory inflammatory foreign body response and transient inflammatory tissue reactions</li><li>• Wound dehiscence</li><li>• Calculi formation in urinary or biliary tracts when prolonged contact with salt solutions such as urine or bile occurs.</li><li>• Like all foreign bodies, FILAPROP may contribute in enhancing an existing infection.</li></ul></div> <div><div>INFORMATION TO PATIENT</div><p>•Keep the wound site dry and clean.</p><p>•Consult the doctor if a suture opens, loosens or breaks.</p><p>•Consult the doctor in case of post operative bleeding or leakage.</p><p>•Consult doctor/ surgeon in case of persistent pain and discomfort, fever, headache, nausea and vomiting.</p></div> <div><div>STERILITY</div><p>FILAPROP™ sutures are sterilized by ethylene oxide (ETO) as indicated on the package. Do not re-sterilize! If re-sterilized, product may degrade and not perform as intended. Do not use if package is opened or damaged. Discard opened remaining unused sutures as per the local area discard procedure !</p></div> <div><div>STORAGE</div><p>Recommended storage conditions: Between 15°C and 30°C, away from moisture and direct heat. Do not use after expiry date.</p></div>
<div><div>FILAPROP™</div><div>NON-ABSORBABLE POLYPROPYLENE SURGICAL SUTURE U.S.P.</div></div> <div><div>DESCRIPTION</div><p>FILAPROP™ suture is a monofilament, synthetic non-absorbable, sterile surgical suture composed of an isotactic crystalline stereoisomer of polypropylene, a synthetic linear polyolefin. The molecular formula of polypropylene is (C<sub>3</sub>H<sub>5</sub>)<sub>n</sub>. FILAPROP™ sutures are either undyed or dyed with phthalocyanine blue.</p><p>FILAPROP™ is available in a range of gauge sizes and lengths attached to standard stainless steel needles of various types and sizes.</p><p>FILAPROP™ complies with the requirements of the United States Pharmacopeia for "Non Absorbable Surgical Suture" and the European Pharmacopeia for "Sterile Synthetic Non Absorbable monofilament Strand." However, it may be slightly oversize in diameter to U.S.P. for some suture sizes. U.S.P except for diameter.</p></div> <div><div>INDICATION</div><p>FILAPROP™ Sutures are indicated for approximation and/or ligation of:</p><ul style="list-style-type: none"><li>- Skin tissue</li><li>- Soft tissues including cardiac</li></ul></div> <div><div>INTENDED USER GROUP</div><p>FILAPROP™ Sutures are intended to be used by Healthcare professionals, general physicians and surgeons who performs skin closures, soft tissue approximation and/or ligation.</p></div> <div><div>INTENDED TARGET GROUP</div><p>FILAPROP™ Suture is intended to be used in both men and women and in all general populations including infants, children, and adults of all age groups.</p></div> <div><div>INTENDED PURPOSE</div><p>FILAPROP™ Sutures are intended to be used for skin closure and soft tissue (internal tissue) approximation and/or ligation.</p></div> <div><div>SELECTION CRITERIA</div><p>The selection of suture for implantation depends on patient condition, surgical experience, surgical technique and wound size.</p><p>The device is used by professional healthcare person only.</p></div> <div><div>PERFORMANCE</div><p>FILAPROP™ suture initially elicits a minimal inflammatory reaction in</p></div>	<div><div>DEVICE LIFETIME</div><p>FILAPROP™ Suture is Long term permanent implantable device.</p></div> <div><div>RISK OF REUSE</div><p>The product is instructed to be used single time only. If used multiple times it may lead to cross contamination / infection.</p></div> <div><div>DISPOSAL OF PRODUCT</div><p>Disposal of unused / used product and packaging should be carried out in accordance with hospital, administrative and / or local government policy.</p></div> <div><div>SSCP SUMMARY</div><p>A summary of safety and clinical performance for FILAPROP™ Suture can be found at the following link (once the EUDAMED will be activated.)</p><p><a href="https://ec.europa.eu/tools/eudamed">https://ec.europa.eu/tools/eudamed</a></p></div> <div><div>"Contact Meril Endo Surgery Pvt. Ltd. and Competent authority of the member state in case of any serious incident occur with the device"</div></div>
<div><div>FILAPROP™</div><div>NON-ABSORBABLE POLYPROPYLENE SURGICAL SUTURE U.S.P.</div></div> <div><div>DESCRIPTION</div><p>FILAPROP™ suture is a monofilament, synthetic non-absorbable, sterile surgical suture composed of an isotactic crystalline stereoisomer of polypropylene, a synthetic linear polyolefin. The molecular formula of polypropylene is (C<sub>3</sub>H<sub>5</sub>)<sub>n</sub>. FILAPROP™ sutures are either undyed or dyed with phthalocyanine blue.</p><p>FILAPROP™ is available in a range of gauge sizes and lengths attached to standard stainless steel needles of various types and sizes.</p><p>FILAPROP™ complies with the requirements of the United States Pharmacopeia for "Non Absorbable Surgical Suture" and the European Pharmacopeia for "Sterile Synthetic Non Absorbable monofilament Strand." However, it may be slightly oversize in diameter to U.S.P. for some suture sizes. U.S.P except for diameter.</p></div> <div><div>INDICATION</div><p>FILAPROP™ Sutures are indicated for approximation and/or ligation of:</p><ul style="list-style-type: none"><li>- Skin tissue</li><li>- Soft tissues including cardiac</li></ul></div> <div><div>INTENDED USER GROUP</div><p>FILAPROP™ Sutures are intended to be used by Healthcare professionals, general physicians and surgeons who performs skin closures, soft tissue approximation and/or ligation.</p></div> <div><div>INTENDED TARGET GROUP</div><p>FILAPROP™ Suture is intended to be used in both men and women and in all general populations including infants, children, and adults of all age groups.</p></div> <div><div>INTENDED PURPOSE</div><p>FILAPROP™ Sutures are intended to be used for skin closure and soft tissue (internal tissue) approximation and/or ligation.</p></div> <div><div>SELECTION CRITERIA</div><p>The selection of suture for implantation depends on patient condition, surgical experience, surgical technique and wound size.</p><p>The device is used by professional healthcare person only.</p></div> <div><div>PERFORMANCE</div><p>FILAPROP™ suture initially elicits a minimal inflammatory reaction in</p></div>	<div><div>SYMBOLS USED ON LABELLING</div><div><div><div>2</div><div></div></div><div>= Do not reuse</div></div><div><div><div>IN</div><div></div></div><div>= Date of Manufacture</div></div><div><div><div></div><div></div></div><div>= Use by (Use until Year &amp; Month)</div></div><div><div><div>STERILE EO</div><div></div></div><div>= Sterilized using Ethylene Oxide</div></div><div><div><div>LOT</div><div></div></div><div>= Batch Code (Number)</div></div><div><div><div></div><div></div></div><div>= Consult Instruction for use</div></div><div><div><div></div><div></div></div><div>= Do not resterilize</div></div><div><div><div>15°C 30°C</div><div></div></div><div>= Temperature Limit (Store between 15°C and 30°C)</div></div><div><div><div></div><div></div></div><div>= Keep away from Sunlight</div></div><div><div><div></div><div></div></div><div>= Keep Dry</div></div><div><div><div></div><div></div></div><div>= Do not use if Package is Damaged</div></div><div><div><div></div><div></div></div><div>= Caution</div></div><div><div><div></div><div></div></div><div>= Manufacturer</div></div><div><div><div></div><div></div></div><div>= Dyed / Monofilament Nonabsorbable</div></div><div><div><div>PP</div><div></div></div><div>= Polypropylene</div></div><div><div><div>EC REP</div><div></div></div><div>= Authorised Representative in the European Community</div></div><div><div><div>CE 0297</div><div></div></div><div>= CE mark and Identification number of Notified Body</div></div><div><div><div>MD</div><div></div></div><div>= Medical Device</div></div><div><div><div></div><div></div></div><div>= Single sterile barrier system with protective packaging outside</div></div><div><div><div>MR</div><div></div></div><div>= MR Safe</div></div><div><div><div>UDI</div><div></div></div><div>= Unique Device Identifier</div></div></div>

ARTWORK No.	REV.
ME/IFU/194	00

PANTONE SHADE
<div></div> PANTONE 278 C
<div></div> BLACK