

Corporate Headquarters, Vapi

Meril Healthcare Pvt. Ltd.
H1-H3, Meril Park
Survey No. 135/2/B & 174/2
Muktanand Marg, Chala
Vapi 396 191. Gujarat. India.
T: +91 260 3052 100
F: +91 260 3052 125
E: orthopedics@merillife.com

International Sales & Marketing

Meril Healthcare Pvt. Ltd.
612, B-Wing, Bonanaza
Sahar Plaza, Andheri East
Mumbai 400 059
Maharashtra. India.
T: +91 22 3935 0700
F: +91 22 4047 9717

India Sales & Marketing

Meril Healthcare Pvt. Ltd.
512, Midas, Sahar Plaza, J.B. Nagar
Andheri East, Mumbai 400 059
Maharashtra. India.
T: +91 22 4047 9797

Meril South America

DOC MED LTDA
1079 – Cep: 04077-003 - Moema
Sao Paulo, Brazil
T/F: +55 11 3624 5935/6

Meril GmbH

Bornheimer Strasse 135 – 137
D – 53119 Bonn, Germany
T/F: +49 228 7100 4000/1

Meril Tibbi Cihazlar

Imalat Ve Ticaret A.S.
6, Mimar Sinan Mah.,
Cavusbasi Cad, Ozde Sok
Aydin Eksi Is Merkezi Kat:1
Cekmekoy/Istanbul, Turkey
T: +90 53 2272 5172



ARMAR™ Catalogue

Titanium Plate | Screw | Instrument

Cardiovascular

Orthopedics

Diagnostics

Endo-Surgery

ENT

Meril

More to Life

More to Life

Meril is a global medical device company dedicated to the design and development of novel, clinically relevant, state-of-the-art and best-in-class devices to alleviate human suffering and improve the quality of life. We span a broad operational canvas ranging from vascular interventional devices to orthopedics, in-vitro diagnostics and endo-surgery.

We share an enduring perseverance to advance healthcare solutions, for more patients to live longer and healthier lives. We have a strong commitment towards R&D and adhere to the best quality standards in manufacturing, scientific communication and distribution known today.

Established in 2006, Meril's launch was in line with the healthcare diversification plan of Bilakhia Group. Located in Vapi, 150km north of Mumbai, and ensconced within the serene hill-sides of Chala town, Meril sits proud in its ultra modern manufacturing facilities.

Where research becomes the basis for discovering new ways of alleviating human suffering and empowering an individual with good health

ri Research
Empowering
Individual




Infrastructure

Situated on 180 acres of land and built over a 300,000 sq.ft. area, Meril's ISO 13485 & cGMP certified ultra-modern manufacturing facility is constructed to satisfy the stringent needs of our cardiovascular, orthopedics, diagnostics and endo-surgery businesses, as well as our R&D endeavors.

Fully integrated manufacturing systems ensure backward integration, flawless man-material movement and complete control over processes to eliminate production errors. All manufacturing and sterilization processes are conducted in-house, in addition to the analytical and microbiological QA/QC tests required to meet the world-class production standards.



Meril's 3,200+ personnel are a strong, experienced team comprising designers, engineers, chemists, microbiologists, regulatory affair experts, R&D scientists, clinical affairs experts, legal, finance, sales and marketing professionals with innovative capabilities. The team is continuously striving towards improving and saving lives every day by providing revolutionary solutions for diagnosis, prevention and treatment.



Global Market Presence

Currently Meril has direct sales & representatives offices in India, Germany, Turkey, USA, China and Brazil and strong distribution model spread over 102 countries across continents.



Meril Orthopedic

Meril is known for an uncompromising attitude when it comes to quality of design. Providing the medical world with zero defect devices is our only motto. Our orthopedic trauma products---Meril Locking Plate System 3.5, Meril Locking Plate System 4.5, Meril Compression Plate System 3.5, Meril Compression Plate System 4.5, Meril Dynamic Hip Fixation System, and Meril Locking Distal Radial System 2.4 are specifically designed to treat fractures, bone tumour and deformities. They boast of high fatigue strength, proper fixation and negligible soft tissue damage.

We have been providing trauma implants for the foot, shoulder, hip, hand, pelvis, and condyles. In order to ensure smooth and trouble-free operating procedures leading medical professionals rely only on our Trauma products. We have carried out performance testing to prove our design for its endurance and durability and the outcome of these tests prove that they are going to create a dynamic image in the history of joint replacement.



Index	
• ARMAR™ 2.4mm Locking Distal Radial Plate System	1
• ARMAR™ 3.5mm Meril Locking Plate System	13
• ARMAR™ Elbow-Clavicle Plate System	33
• ARMAR™ 4.5mm Meril Locking Plate System	47
• MBOSS™ 4.5mm Cannulated Screw	63
• MBOSS™ 6.5mm Cannulated Screw	69
• ARMAR™ Proximal Humeral Locking Plate System	77
• ARMAR™ Proximal Lateral Tibia Locking Plate System	81
• ARMAR™ Distal Femoral Locking Plate System	87

ARMAR™ 2.4 Locking Distal Radial Plate System

LPS Distal Radial Plate 2.4, Volar | Material : TiCP | Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PT5303043L	5	3	43
MT-PT5303043R	5	3	43
MT-PT5303061L	5	5	61
MT-PT5303061R	5	5	61



LPS Distal Radial Plate 2.4, Volar Buttress | Material : TiCP | Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PT5304042L	5	3	42
MT-PT5304042R	5	3	42
MT-PT5304057L	5	4	52
MT-PT5304057R	5	4	52



LPS Distal Radial Plate 2.4, Extra-articular | Material : TiCP | Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PT5303048L	5	3	48
MT-PT5305066L	5	5	66
MT-PT5303047L	4	3	47
MT-PT5305065L	4	5	65
MT-PT5303048R	5	3	48
MT-PT5305066R	5	5	66
MT-PT5303047R	4	3	47
MT-PT5305065R	4	5	65



LPS T-Distal Radial Plate 2.4 | Material : TiCP | Thickness : 1.6mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PT4203040	3	3	40
MT-PT4204049	3	4	49



ARMAR™ 2.4 Locking Distal Radial Plate System

LPS Distal Radial Plate 2.4, Straight | Material : TiCP | Thickness : 1.6mm

PRODUCT CODE	SHAFT HOLES	LENGTH
MT-PT4205048	5	48
MT-PT4206057	6	57



LPS L-Distal Radial Plate 2.4 Right Angled | Material : TiCP | Thickness : 1.6mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PT4202040L	2	3	40
MT-PT4202049L	2	4	49
MT-PT4203040L	3	3	40
MT-PT4203049L	3	4	49
MT-PT4202040R	2	3	40
MT-PT4202049R	2	4	49
MT-PT4203040R	3	3	40
MT-PT4203049R	3	4	49



LPS L-Distal Radial Plate 2.4, Oblique angled | Material : TiCP | Thickness : 1.6mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PT4203043L	3	3	43
MT-PT4204052L	3	4	52
MT-PT4203043R	3	3	43
MT-PT4204052R	3	4	52



LPS Distal Radial Plate 2.4, Long | Material : TiCP | Thickness : 1.6mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PT5308100	4	8	100
MT-PT5310120	4	10	120
MT-PT5312140	4	12	140



ARMAR™ 2.4 Locking Distal Radial Plate System

2.4mm LPS Volar Rim Distal Radial Plate | Material : TiCP | Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PT7106057R	6	5	57
MT-PT7106057L	6	5	57
MT-PT7107057R	7	5	57
MT-PT7107057L	7	5	57



2.4mm LPS Double-Column Distal Radial Plate, Volar 19.5 mm

Material : TiCP | Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PT6802042R	6	2	42
MT-PT6803051R	6	3	51
MT-PT6804063R	6	4	63
MT-PT6805072R	6	5	72
MT-PT6802042L	6	2	42
MT-PT6803051L	6	3	51
MT-PT6804063L	6	4	63
MT-PT6805072L	6	5	72



2.4mm LPS Double-Column Distal Radial Plate, Volar 22.5 mm

Material : TiCP | Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PT6902045R	6	2	45
MT-PT6903054R	6	3	54
MT-PT6904066R	6	4	66
MT-PT6905075R	6	5	75
MT-PT6902045L	6	2	45
MT-PT6903054L	6	3	54
MT-PT6904066L	6	4	66
MT-PT6905075L	6	5	75



ARMAR™ 2.4 Locking Distal Radial Plate System

2.4mm LPS Double-Column Distal Radial Plate, Volar 25.5 mm

Material : TiCP | Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PT7002047R	7	2	47
MT-PT7003056R	7	3	56
MT-PT7004068R	7	4	68
MT-PT7005077R	7	5	77
MT-PT7002047L	7	2	47
MT-PT7003056L	7	3	56
MT-PT7004068L	7	4	68
MT-PT7005077L	7	5	77



MBOSS™ 2.4mm Screw

Ø 2.4mm Locking Screw T8, Self Tapping | Material : TAN
Head Diameter : 3.6mm | Outer Diameter : 2.4mm | Core Diameter : 1.9mm | Pitch : 0.6mm

PRODUCT CODE

LENGTH

MT-ST0124006	6
MT-ST0124007	7
MT-ST0124008	8
MT-ST0124009	9
MT-ST0124010	10
MT-ST0124011	11
MT-ST0124012	12
MT-ST0124013	13
MT-ST0124014	14
MT-ST0124016	16
MT-ST0124018	18
MT-ST0124020	20
MT-ST0124022	22
MT-ST0124024	24
MT-ST0124026	26
MT-ST0124028	28
MT-ST0124030	30



Ø 2.4mm Cortex Screw T8, Self Tapping | Material : TAN
Head Diameter : 4.0mm | Outer Diameter : 2.4mm | Core Diameter : 1.7mm | Pitch : 1.0mm

PRODUCT CODE

LENGTH

MT-ST0224006	6
MT-ST0224007	7
MT-ST0224008	8
MT-ST0224009	9
MT-ST0224010	10
MT-ST0224012	12
MT-ST0224014	14
MT-ST0224016	16
MT-ST0224018	18
MT-ST0224020	20
MT-ST0224022	22
MT-ST0224024	24
MT-ST0224026	26
MT-ST0224028	28
MT-ST0224030	30
MT-ST0224032	32
MT-ST0224034	34
MT-ST0224036	36
MT-ST0224038	38
MT-ST0224040	40



MBOSS™ 2.4mm Screw

Ø 2.7mm Cortex Screw T8, Self Tapping | Material : TAN
Head Diameter : 5.0mm | Outer Diameter : 2.7mm | Core Diameter : 1.9mm | Pitch : 1.0mm

PRODUCT CODE

LENGTH

MT-ST0527008	8
MT-ST0527010	10
MT-ST0527012	12
MT-ST0527014	14
MT-ST0527016	16
MT-ST0527018	18
MT-ST0527020	20
MT-ST0527022	22
MT-ST0527024	24
MT-ST0527026	26
MT-ST0527028	28
MT-ST0527030	30
MT-ST0527032	32
MT-ST0527034	34
MT-ST0527036	36
MT-ST0527038	38
MT-ST0527040	40
MT-ST0527045	45
MT-ST0527050	50
MT-ST0527055	55
MT-ST0527060	60



MBOSS™ 2.4mm Screw

2.4mm Headless Compression Screw Short Thread

Material : TAN | Head Diameter : 3.1 | Outer Diameter : 2.4 | Core Diameter : 1.9 | Pitch : 1.0

PRODUCT CODE

LENGTH

MT-ST0324009	09
MT-ST0324010	10
MT-ST0324012	12
MT-ST0324014	14
MT-ST0324016	16
MT-ST0324018	18
MT-ST0324020	20
MT-ST0324022	22
MT-ST0324024	24
MT-ST0324026	26
MT-ST0324028	28
MT-ST0324030	30
MT-ST0324032	32
MT-ST0324034	34
MT-ST0324036	36
MT-ST0324038	38
MT-ST0324040	40



2.4mm Headless Compression Screw Long Thread

Material : TAN | Head Diameter : 3.1 | Outer Diameter : 2.4 | Core Diameter : 1.9 | Pitch : 1.0


PRODUCT CODE

LENGTH

MT-ST0424016	16
MT-ST0424018	18
MT-ST0424020	20
MT-ST0424022	22
MT-ST0424024	24
MT-ST0424026	26
MT-ST0424028	28
MT-ST0424030	30



ARMAR™ 2.4mm Locking Plate System Instrument Set



ARMAR™ 2.4mm Locking Plate System Instrument Set

Part No:MT-INDR24001

Description

Drill bit Ø1.8mm, Length 110 2-flute, for Quick Coupling



Part No:MT-INDR24002

Description


Drill Bit Ø2.0mm, Length 110 2-flute, for Quick Coupling



Part No:MT-INDR24007

Description

Universal Drill Guide 2.4



Part No:MT-INDR24008

Description

Universal Drill Guide 2.7



Part No:MT-INDR24005

Description


LPS Drill Sleeve 2.4, with Scale up to 30mm, for Drill bits Ø1.8mm



ARMAR™ 2.4mm Locking Plate System Instrument Set

Part No:MT-INDR24006	Description
	LPS Drill Sleeve 2.7 with Scale up to 30mm, for Drill bits Ø2.0mm
	
Part No:MT-INDR24015	Description
	Screwdriver shaft, T8, self-holding
	
Part No:MT-INDR24013	Description
	Manual 0.8Nm Torque Limiter
	
Part No:MT-INDR24020	Description
	Quick Change Handle Large with Quick Coupling, Length 138
	
Part No:MT-INDR24018	Description
	Depth Gauge for Screws 2.0/2.4mm, measuring range up to 40mm
	
Part No:MT-INDR24021	Description
	Grasping sets of 2.4/2.7 locking screw, Screw Driver T8
	

ARMAR™ 2.4mm Locking Plate System Instrument Set

Part No:MT-INDR24016	Description
	Screw Driver Stardrive T8 Self Holding
	
Part No:MT-INDR24028	Description
	Ø1.25mm Guide Wire Length 150mm
	
Part No:MT-INDR24027	Description
	Instrument Case with Lid-DR24
	

ARMAR™ 3.5mm Meril Locking Plate System

LPS-Small Compression Plate | Material : TiCP | Thickness : 3.3mm

PRODUCT CODE	HOLES	LENGTH
MT-PT0804059	4	59
MT-PT0805072	5	72
MT-PT0806085	6	85
MT-PT0807098	7	98
MT-PT0808111	8	111
MT-PT0809124	9	124
MT-PT0810137	10	137
MT-PT0811150	11	150
MT-PT0812163	12	163
MT-PT0813176	13	176
MT-PT0814189	14	189
MT-PT0815202	15	202
MT-PT0816215	16	215
MT-PT0818241	18	241
MT-PT0820267	20	267
MT-PT0822293	22	293



LPS-Small Metaphyseal Plate | Material : TiCP | Thickness : 3.3mm

PRODUCT CODE	HOLES	LENGTH
MT-PT1606086	6	86
MT-PT1607099	7	99
MT-PT1608112	8	112
MT-PT1609125	9	125
MT-PT1610138	10	138
MT-PT1611151	11	151
MT-PT1612164	12	164
MT-PT1614190	14	190
MT-PT1616216	16	216
MT-PT1618242	18	242



ARMAR™ 3.5mm Meril Locking Plate System

LPS-Reconstruction Round Hole Plate 3.5 | Material : TiCP | Thickness : 3.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PT0705058	5	58
MT-PT0706070	6	70
MT-PT0707082	7	82
MT-PT0708094	8	94
MT-PT0709106	9	106
MT-PT0710118	10	118
MT-PT0712142	12	142
MT-PT0714166	14	166



LPS-Reconstruction Combined Hole Plate 3.5 | Material : TiCP | Thickness : 3.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PT1804056	4	56
MT-PT1805070	5	70
MT-PT1806084	6	84
MT-PT1807098	7	98
MT-PT1808112	8	112
MT-PT1809126	9	126
MT-PT1810140	10	140
MT-PT1811154	11	154
MT-PT1812168	12	168
MT-PT1814196	14	196
MT-PT1816224	16	224
MT-PT1818252	18	252
MT-PT1820280	20	280
MT-PT1822308	22	308



ARMAR™ 3.5mm Meril Locking Plate System


LPS-T Plate Right Angled 3.5 | Material : TiCP | Thickness : 1.8mm

PRODUCT CODE	HOLES		LENGTH
	HEAD	SHAFT	
MT-PT1903050R	3	3	46
MT-PT1904056R	3	4	56.5
MT-PT1905067R	3	5	67
MT-PT1904056R	4	4	57
MT-PT1906078R	4	6	78



LPS Cloverleaf Plate | Material : TiCP | Thickness : 2.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PT1303088	3	92
MT-PT1304104	4	107
MT-PT1305120	5	121
MT-PT1306136	6	136



LPS T Plate Oblique Angled 3.5 | Material : TiCP | Thickness : 1.8mm

PRODUCT CODE	HOLES		LENGTH
	HEAD	SHAFT	
MT-PT2003052L	3	3	52
MT-PT2004063L	3	4	63
MT-PT2005074L	3	5	74
MT-PT2006085L	3	6	85
MT-PT2007096L	3	7	96
MT-PT2008107L	3	8	107
MT-PT2003052R	3	3	52
MT-PT2004063R	3	4	63
MT-PT2005074R	3	5	74
MT-PT2006085R	3	6	85
MT-PT2007096R	3	7	96
MT-PT2008107R	3	8	107



ARMAR™ 3.5mm Meril Locking Plate System

LPS Proximal Humerus Plate Short 3.5 | Material : TiCP | Thickness : 3.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PT5103090	3	90
MT-PT5105114	5	114



LPS Proximal Humerus Plate Long 3.5 | Material : TiCP | Thickness : 4.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PT5205142	5	142
MT-PT5206160	6	160
MT-PT5208196	8	196
MT-PT5210232	10	232
MT-PT5212268	12	268



LPS Proximal Tibia Plate 3.5 | Material : TiCP | Thickness : 4.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PT1704081L	4	76
MT-PT1706107L	6	102
MT-PT1708133L	8	128
MT-PT1710159L	10	154
MT-PT1712185L	12	180
MT-PT1714211L	14	206
MT-PT1716237L	16	232
MT-PT1704081R	4	76
MT-PT1706107R	6	102
MT-PT1708133R	8	128
MT-PT1710159R	10	154
MT-PT1712185R	12	180
MT-PT1714211R	14	206
MT-PT1716237R	16	232



ARMAR™ 3.5mm Meril Locking Plate System

LPS Medial Proximal Tibial Plate 3.5 | Material : TiCP | Thickness : 4.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PT1404093L	4	93
MT-PT1406119L	6	119
MT-PT1408145L	8	145
MT-PT1410171L	10	171
MT-PT1412197L	12	197
MT-PT1414223L	14	223
MT-PT1416249L	16	249
MT-PT1418275L	18	275
MT-PT1420301L	20	301
MT-PT1404093R	4	93
MT-PT1406119R	6	119
MT-PT1408145R	8	145
MT-PT1410171R	10	171
MT-PT1412197R	12	197
MT-PT1414223R	14	223
MT-PT1416249R	16	249
MT-PT1418275R	18	275
MT-PT1420301R	20	301



LPS Medial Distal Tibia Low Bend Plate 3.5 | Material : TiCP | Thickness : 4.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PT4704109L	4	109
MT-PT4706135L	6	135
MT-PT4708161L	8	161
MT-PT4710187L	10	187
MT-PT4712213L	12	213
MT-PT4714239L	14	239
MT-PT4704109R	4	109
MT-PT4706135R	6	135
MT-PT4708161R	8	161
MT-PT4710187R	10	187
MT-PT4712213R	12	213
MT-PT4714239R	14	239



ARMAR™ 3.5mm Meril Locking Plate System

LPS One Third Tubular Plate 3.5 | Material : TiCP | Thickness : 1.2mm

PRODUCT CODE	HOLES	LENGTH
MT-PT5002028	2	28
MT-PT5003040	3	40
MT-PT5004052	4	52
MT-PT5005064	5	64
MT-PT5006076	6	76
MT-PT5007088	7	88
MT-PT5008100	8	100
MT-PT5009112	9	112
MT-PT5010124	10	124
MT-PT5012148	12	148




LPS Proximal Humerus Locking Plate 3.5 | Material : TiCP | Thickness : 2.5mm

PRODUCT CODE	HOLES	LENGTH
MT-PT4905085	5	85
MT-PT4908121	8	121



LPS Anterior-Lateral Distal Tibial Plate | Material : TiCP | Thickness : 4.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PT1205080L	5	80
MT-PT1207106L	7	106
MT-PT1209132L	9	132
MT-PT1211158L	11	158
MT-PT1213184L	13	184
MT-PT1215210L	15	210
MT-PT1217236L	17	236
MT-PT1219262L	19	262
MT-PT1221288L	21	288
MT-PT1205080R	5	80
MT-PT1207106R	7	106
MT-PT1209132R	9	132
MT-PT1211158R	11	158
MT-PT1213184R	13	184
MT-PT1215210R	15	210
MT-PT1217236R	17	236
MT-PT1219262R	19	262
MT-PT1221288R	21	288



ARMAR™ 3.5mm Meril Locking Plate System

3.5mm LPS Calcaneal Plate

Material : TiCP | Thickness : 2.5 mm

PRODUCT CODE	TYPE	LENGTH
MT-PT2115064L	Mini	64
MT-PT2115069L	Small	69
MT-PT2115076L	Long	76
MT-PT2115064R	Mini	64
MT-PT2115069R	Small	69
MT-PT2115076R	Long	76



3.5mm LPS Periarticular Proximal Lateral Tibia Plate

Material : TiCP | Thickness : 3.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PT7306112L	6	112
MT-PT7308138L	8	138
MT-PT7310164L	10	164
MT-PT7312190L	12	190
MT-PT7314216L	14	216
MT-PT7316242L	16	242
MT-PT7306112R	6	112
MT-PT7308138R	8	138
MT-PT7310164R	10	164
MT-PT7312190R	12	190
MT-PT7314216R	14	216
MT-PT7316242R	16	242



3.5mm LPS Periarticular Proximal Lateral Humeral Plate

Material : TiCP | Thickness : 3.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PT7404078L	4	78
MT-PT7406104L	6	104
MT-PT7408130L	8	130
MT-PT7412182L	12	182
MT-PT7416234L	16	234
MT-PT7404078R	4	78
MT-PT7406104R	6	104
MT-PT7408130R	8	130
MT-PT7412182R	12	182
MT-PT7416234R	16	234



ARMAR™ 3.5mm Meril Locking Plate System

3.5mm LPS Proximal Tibia Plate, Posteromedial

Material : TiCP | Thickness : 3.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PT7501069	1	69
MT-PT7502079	2	79
MT-PT7504105	4	105
MT-PT7506131	6	131
MT-PT7508157	8	157
MT-PT7510183	10	183



2.7mm/3.5mm LPS Lateral Distal Fibula Plate

Material : TiCP | Thickness : 3.0 mm

PRODUCT CODE	HOLES	LENGTH
MT-PT6704080R	4	80
MT-PT6706106R	6	106
MT-PT6708132R	8	132
MT-PT6710158R	10	158
MT-PT6704080L	4	80
MT-PT6706106L	6	106
MT-PT6708132L	8	132
MT-PT6710158L	10	158




MBOSS™ 3.5mm Screw

Ø 3.5mm Locking Screw, Self Tapping | Material : TAN
Head Diameter : 5.0mm | Outer Diameter : 3.5mm | Core Diameter : 2.9mm | Pitch : 0.8mm

PRODUCT CODE

LENGTH

MT-ST1235010	10
MT-ST1235012	12
MT-ST1235014	14
MT-ST1235016	16
MT-ST1235018	18
MT-ST1235020	20
MT-ST1235022	22
MT-ST1235024	24
MT-ST1235026	26
MT-ST1235028	28
MT-ST1235030	30
MT-ST1235032	32
MT-ST1235034	34
MT-ST1235036	36
MT-ST1235038	38
MT-ST1235040	40
MT-ST1235042	42
MT-ST1235044	44
MT-ST1235046	46
MT-ST1235048	48
MT-ST1235050	50
MT-ST1235052	52
MT-ST1235054	54
MT-ST1235056	56
MT-ST1235058	58
MT-ST1235060	60
MT-ST1235062	62
MT-ST1235065	65
MT-ST1235070	70
MT-ST1235075	75
MT-ST1235080	80
MT-ST1235085	85
MT-ST1235090	90
MT-ST1235095	95



MBOSS™ 3.5mm Screw

Ø 3.5mm Cortex Screw, Self Tapping | Material : TiCP
Head Diameter : 6.0mm | Outer Diameter : 3.5mm | Core Diameter : 2.4mm | Pitch : 1.25mm

PRODUCT CODE

LENGTH

MT-ST1135010	10
MT-ST1135012	12
MT-ST1135014	14
MT-ST1135016	16
MT-ST1135018	18
MT-ST1135020	20
MT-ST1135022	22
MT-ST1135024	24
MT-ST1135026	26
MT-ST1135028	28
MT-ST1135030	30
MT-ST1135032	32
MT-ST1135034	34
MT-ST1135036	36
MT-ST1135038	38
MT-ST1135040	40
MT-ST1135042	42
MT-ST1135044	44
MT-ST1135046	46
MT-ST1135048	48
MT-ST1135050	50
MT-ST1135055	55
MT-ST1135060	60
MT-ST1135065	65
MT-ST1135070	70
MT-ST1135075	75
MT-ST1135080	80
MT-ST1135085	85
MT-ST1135090	90



MBOSS™ 3.5mm Screw

Ø 4.0mm Cancellous Screw, Short Thread | Material : TiCP
Head Diameter : 6.0mm | Outer Diameter : 4.0mm | Core Diameter : 2.0mm | Pitch : 1.75mm

PRODUCT CODE

LENGTH

MT-ST1940010	10
MT-ST1940012	12
MT-ST1940014	14
MT-ST1940016	16
MT-ST1940018	18
MT-ST1940020	20
MT-ST1940022	22
MT-ST1940024	24
MT-ST1940026	26
MT-ST1940028	28
MT-ST1940030	30
MT-ST1940035	35
MT-ST1940040	40
MT-ST1940045	45
MT-ST1940050	50
MT-ST1940055	55
MT-ST1940060	60
MT-ST1940065	65
MT-ST1940070	70
MT-ST1940075	75
MT-ST1940080	80



MBOSS™ 3.5mm Screw

Ø 4.0mm Cancellous Screw, Full Thread Screw | Material : TiCP
Head Diameter : 6.0mm | Outer Diameter : 4.0mm | Core Diameter : 2.0mm | Pitch : 1.75mm

PRODUCT CODE

LENGTH

MT-ST1840010	10
MT-ST1840012	12
MT-ST1840014	14
MT-ST1840016	16
MT-ST1840018	18
MT-ST1840020	20
MT-ST1840022	22
MT-ST1840024	24
MT-ST1840026	26
MT-ST1840028	28
MT-ST1840030	30
MT-ST1840040	35
MT-ST1840035	40
MT-ST1840045	45
MT-ST1840050	50
MT-ST1840055	55
MT-ST1840060	60
MT-ST1840065	65
MT-ST1840070	70
MT-ST1840075	75
MT-ST1840080	80



MBOSS™ 3.5mm Screw

Ø 3.5mm Locking Cancellous Screw, Self Tapping | Material : TAN
Head Diameter : 5.0 mm | Outer Diameter : 3.5 mm | Core Diameter : 2.0 mm | Pitch : 1.75mm

PRODUCT CODE

LENGTH

MT-ST8435030	30
MT-ST8435035	35
MT-ST8435040	40
MT-ST8435045	45
MT-ST8435050	50
MT-ST8435055	55
MT-ST8435060	60
MT-ST8435065	65
MT-ST8435070	70
MT-ST8435075	75
MT-ST8435080	80
MT-ST8435085	85
MT-ST8435090	90



ARMAR™ 3.5mm Locking Plate System Instrument Set



ARMAR™ 3.5mm Locking Plate System Instrument Set

Part No:MT-INSLPS002	Description
	Drill Bit Ø2.5mm, Length 125mm, 2 flute, for Quick coupling
	
Part No:MT-INSLPS003	Description
	Drill Bit Ø2.8mm, Length 147mm, 2-flute, for Quick coupling
	
Part No:MT-INSLPS004	Description
	Drill Bit Ø3.5mm, Length 125mm, 2 flute, for Quick coupling
	
Part No:MT-INSLPS044	Description
	Drill Bit Ø2.8mm, Length 200mm, 2-flute, for Quick coupling
	
Part No:MT-INSLPS006	Description
	Tap for Cortex screw Ø3.5mm Length 115mm for Quick Coupling
	
Part No:MT-INSLPS020	Description
	LPS Drill Sleeve for Drill Bits Ø2.8mm, Length 90mm
	

ARMAR™ 3.5mm Locking Plate System Instrument Set

Part No:MT-INSLPS012

Description

LPS Drill Sleeve for Drill Bits Ø2.8mm, Length 65mm



Part No:MT-INSLPS015

Description

Universal Drill Guide 2.5mm & 3.5mm



Part No:MT-INSLPS011

Description

Double drill Guide 3.5mm & 2.5 mm



Part No:MT-INSLPS023

Description

2.5 Screw Driver shaft with quick coupling for Torque Limiter



Part No:MT-INSLPS024

Description

Screwdriver, Hexagonal, small, A/F 2.5mm, with Groove & Holding Sleeve



ARMAR™ 3.5mm Locking Plate System Instrument Set

Part No:MT-INSLPS022

Description

Manual 1.5Nm Torque Limiter



Part No:MT-INSLPS026

Description

Wrench-SLPS



Part No:MT-INSLPS032

Description

Quick Change Handle for SLPS



Part No:MT-INSLPS031

Description

T-Handle with Quick coupling SLPS



ARMAR™ 3.5mm Locking Plate System Instrument Set

Part No:MT-INSLPS028

Description


Depth gauge for Screw 2.7/3.5mm measuring range up to 90mm



Part No:MT-INSLPS040

Description

Instrument Case with Lid - SLPS



ARMAR™ Elbow-Clavicle Plate System

2.7mm/3.5mm Superior Anterior Clavicle Plate

Material : TiCP | Thickness : 3 mm

PRODUCT CODE	HOLES	LENGTH
MT-PT1006094L	6	92
MT-PT1007110L	7	106
MT-PT1008120L	8	120
MT-PT1006094R	6	92
MT-PT1007110R	7	106
MT-PT1008120R	8	120



2.7mm/3.5mm Superior Anterior Clavicle Plate, with Lateral Extension

Material : TiCP | Thickness : 3 mm

PRODUCT CODE	HOLES	LENGTH
MT-PT1103069L	3	66
MT-PT1104081L	4	80
MT-PT1105094L	5	94
MT-PT1106108L	6	108
MT-PT1107122L	7	122
MT-PT1108136L	8	136
MT-PT1103069R	3	66
MT-PT1104081R	4	80
MT-PT1105094R	5	94
MT-PT1106108R	6	108
MT-PT1107122R	7	122
MT-PT1108136R	8	136



2.7mm/3.5mm Posterolateral Distal Humerus Plate

Material : TiCP | Thickness : 2.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PT0203055L	3	65
MT-PT0205090L	5	91
MT-PT0207116L	7	117
MT-PT0209143L	9	143
MT-PT0214208L	14	208
MT-PT0203055R	3	65
MT-PT0205090R	5	91
MT-PT0207116R	7	117
MT-PT0209143R	9	143
MT-PT0214208R	14	208



ARMAR™ Elbow-Clavicle Plate System

2.7mm/3.5mm Posterolateral Distal Humerus Plate, with Lateral Support

Material : TiCP | Thickness : 2.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PT0303055L	3	65
MT-PT0305090L	5	91
MT-PT0307116L	7	117
MT-PT0309143L	9	143
MT-PT0314208L	14	208
MT-PT0303055R	3	65
MT-PT0305090R	5	91
MT-PT0307116R	7	117
MT-PT0309143R	9	143
MT-PT0314208R	14	208



2.7mm/3.5mm Medial Distal Humerus Plate

Material : TiCP | Thickness : 2.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PT0103058L	3	58
MT-PT0105083L	5	84
MT-PT0107110L	7	110
MT-PT0109149L	9	136
MT-PT0103058R	3	58
MT-PT0105083R	5	84
MT-PT0107110R	7	110
MT-PT0109149R	9	136



2.7mm/3.5mm LPS Medial Distal Humerus Metaphyseal Plate

Material : TiCP | Thickness : 2.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PT1507080	7	80
MT-PT1509106	9	106
MT-PT1511132	11	132
MT-PT1513158	13	158
MT-PT1515184	15	184



ARMAR™ Elbow-Clavicle Plate System

LPS Distal Humerus Extra Articular Plate 3.5

Material : TiCP | Thickness : 2.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PT4104122L	4	122
MT-PT4106158L	6	158
MT-PT4108194L	8	194
MT-PT4110230L	10	230
MT-PT4112266L	12	266
MT-PT4114302L	14	302
MT-PT4104122R	4	122
MT-PT4106158R	6	158
MT-PT4108194R	8	194
MT-PT4110230R	10	230
MT-PT4112266R	12	266
MT-PT4114302R	14	302



LPS Olecranon Plate 3.5

Material : TiCP | Thickness : 3.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PT0902086L	2	86
MT-PT0904112L	4	112
MT-PT0906138L	6	138
MT-PT0908164L	8	164
MT-PT0910190L	10	190
MT-PT0912216L	12	216
MT-PT0902086R	2	86
MT-PT0904112R	4	112
MT-PT0906138R	6	138
MT-PT0908164R	8	164
MT-PT0910190R	10	190
MT-PT0912216R	12	216



ARMAR™ Elbow-Clavicle Plate System

LPS Clavicle Hook Plate 3.5 | Material : TiCP | Thickness : 4.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PT4404012L	4	12mm Depth
MT-PT4405012L	5	
MT-PT4406012L	6	
MT-PT4407012L	7	15mm Depth
MT-PT4404015L	4	
MT-PT4405015L	5	
MT-PT4406015L	6	18mm Depth
MT-PT4407015L	7	
MT-PT4404018L	4	
MT-PT4405018L	5	12mm Depth
MT-PT4406018L	6	
MT-PT4407018L	7	
MT-PT4404012R	4	15mm Depth
MT-PT4405012R	5	
MT-PT4406012R	6	
MT-PT4407012R	7	18mm Depth
MT-PT4404015R	4	
MT-PT4405015R	5	
MT-PT4406015R	6	12mm Depth
MT-PT4407015R	7	
MT-PT4404018R	4	
MT-PT4405018R	5	15mm Depth
MT-PT4406018R	6	
MT-PT4407018R	7	




MBOSS™ 3.5mm Screw

Ø 3.5mm Locking Screw, Self Tapping | Material : TAN
Head Diameter : 5.0mm | Outer Diameter : 3.5mm | Core Diameter : 2.9mm | Pitch : 0.8mm

PRODUCT CODE

LENGTH

MT-ST1235010	10
MT-ST1235012	12
MT-ST1235014	14
MT-ST1235016	16
MT-ST1235018	18
MT-ST1235020	20
MT-ST1235022	22
MT-ST1235024	24
MT-ST1235026	26
MT-ST1235028	28
MT-ST1235030	30
MT-ST1235032	32
MT-ST1235034	34
MT-ST1235036	36
MT-ST1235038	38
MT-ST1235040	40
MT-ST1235042	42
MT-ST1235044	44
MT-ST1235046	46
MT-ST1235048	48
MT-ST1235050	50
MT-ST1235052	52
MT-ST1235054	54
MT-ST1235056	56
MT-ST1235058	58
MT-ST1235060	60
MT-ST1235062	62
MT-ST1235065	65
MT-ST1235070	70
MT-ST1235075	75
MT-ST1235080	80
MT-ST1235085	85
MT-ST1235090	90



MBOSS™ 3.5mm Screw

Ø 3.5mm Cortex Screw, Self Tapping | Material : TiCP
Head Diameter : 6.0 | Outer Diameter : 3.5mm | Core Diameter : 2.4mm | Pitch : 1.25mm

PRODUCT CODE

LENGTH

MT-ST1135010	10
MT-ST1135012	12
MT-ST1135014	14
MT-ST1135016	16
MT-ST1135018	18
MT-ST1135020	20
MT-ST1135022	22
MT-ST1135024	24
MT-ST1135026	26
MT-ST1135028	28
MT-ST1135030	30
MT-ST1135032	32
MT-ST1135034	34
MT-ST1135036	36
MT-ST1135038	38
MT-ST1135040	40
MT-ST1135042	42
MT-ST1135044	44
MT-ST1135046	46
MT-ST1135048	48
MT-ST1135050	50
MT-ST1135055	55
MT-ST1135060	60
MT-ST1135065	65
MT-ST1135070	70
MT-ST1135075	75
MT-ST1135080	80
MT-ST1135085	85
MT-ST1135090	90



MBOSS™ 3.5mm Screw

Ø 4.0mm Cancellous Screw, Full Thread | Material : TiCP | Diameter : 6.0mm
Outer Diameter : 4.0mm | Core Diameter : 2.0mm | Pitch : 1.75mm

PRODUCT CODE	LENGTH
MT-ST1840010	10
MT-ST1840012	12
MT-ST1840014	14
MT-ST1840016	16
MT-ST1840018	18
MT-ST1840020	20
MT-ST1840022	22
MT-ST1840024	24
MT-ST1840026	26
MT-ST1840028	28
MT-ST1840030	30
MT-ST1840035	35
MT-ST1840040	40
MT-ST1840045	45
MT-ST1840050	50
MT-ST1840055	55
MT-ST1840060	60
MT-ST1840065	65
MT-ST1840070	70
MT-ST1840075	75
MT-ST1840080	80



MBOSS™ 3.5mm Screw

Ø 4.0mm Cancellous Screw, Short Thread | Material : TiCP | Diameter : 6.0mm
Outer Diameter : 4.0mm | Core Diameter : 2.0mm | Pitch : 1.75mm

PRODUCT CODE	LENGTH
MT-ST1940010	10
MT-ST1940012	12
MT-ST1940014	14
MT-ST1940016	16
MT-ST1940018	18
MT-ST1940020	20
MT-ST1940022	22
MT-ST1940024	24
MT-ST1940026	26
MT-ST1940028	28
MT-ST1940030	30
MT-ST1940035	35
MT-ST1940040	40
MT-ST1940045	45
MT-ST1940050	50
MT-ST1940055	55
MT-ST1940060	60
MT-ST1940065	65
MT-ST1940070	70
MT-ST1940075	75
MT-ST1940080	80



MBOSS™ 2.7mm Screw

Ø 2.7mm Locking Screw, T8 Self Tapping | Material : TAN
Head Diameter : 3.5 | Outer Diameter : 2.7mm | Core Diameter : 2.1mm | Pitch : 0.6mm

PRODUCT CODE

LENGTH

MT-ST0727012	12
MT-ST0727014	14
MT-ST0727016	16
MT-ST0727018	18
MT-ST0727020	20
MT-ST0727022	22
MT-ST0727024	24
MT-ST0727026	26
MT-ST0727028	28
MT-ST0727030	30
MT-ST0727032	32
MT-ST0727034	34
MT-ST0727036	36
MT-ST0727038	38
MT-ST0727040	40
MT-ST0727042	42
MT-ST0727044	44
MT-ST0727046	46
MT-ST0727048	48
MT-ST0727050	50
MT-ST0727055	55
MT-ST0727060	60



ARMAR™ Elbow and Clavicle System Instrument Set



ARMAR™ Elbow and Clavicle System Instrument Set

Part No:MT-INDR24002

Description

Drill Bit Ø2.0mm, Length 110mm, 2-flute, for Quick Coupling



Part No:MT-INDR24015

Description

Screwdriver shaft, T8, self-holding



Part No:MT-INDR24013

Description

Manual 0.8Nm Torque Limiter



Part No:MT-INDR24020

Description

Quick Change Handle Large with Quick Coupling, Length 138



Part No:MT-INDR24021

Description

Grasping sets of 2.4/2.7 locking screw, Screw Driver T8



Part No:MT-INDR24006

Description

LPS Drill Sleeve 2.7 with Scale up to 30mm, for Drill bits Ø2.0mm



ARMAR™ Elbow and Clavicle System Instrument Set

Part No:MT-INDR24016

Description

Screw Driver Stardrive T8 Self Holding



Part No:MT-INDR24019

Description

Depth Gauge for Screws Ø2.7 to 4.0 mm, measuring range up to 60 mm



Part No:MT-INELCV001

Description

Instrument Case with Lid



ARMAR™ 4.5mm Meril Locking Plate System

LPS Compression Plate Narrow 4.5 | Material : TiCP | Thickness : 4.2mm

PRODUCT CODE	HOLES	LENGTH
MT-PT2902044	2	44
MT-PT2903062	3	62
MT-PT2904080	4	80
MT-PT2905098	5	98
MT-PT2906116	6	116
MT-PT2907134	7	134
MT-PT2908152	8	152
MT-PT2909170	9	170
MT-PT2910188	10	188
MT-PT2911206	11	206
MT-PT2912224	12	224
MT-PT2914260	14	260
MT-PT2916296	16	296



LPS Compression Plate Broad 4.5 | Material : TiCP | Thickness : 5.2mm

PRODUCT CODE	HOLES	LENGTH
MT-PT2606116	6	116
MT-PT2607134	7	134
MT-PT2608152	8	152
MT-PT2609170	9	170
MT-PT2610188	10	188
MT-PT2611206	11	206
MT-PT2612224	12	224
MT-PT2613242	13	242
MT-PT2614260	14	260
MT-PT2615278	15	278
MT-PT2616296	16	296
MT-PT2617314	17	314
MT-PT2618332	18	332



ARMAR™ 4.5mm Meril Locking Plate System

LPS Distal Femur Plate 4.5 | Material : TAN | Thickness : 6.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PT4605156L	5	156
MT-PT4607196L	7	196
MT-PT4609236L	9	236
MT-PT4611276L	11	276
MT-PT4613316L	13	316
MT-PT4605156R	5	156
MT-PT4607196R	7	196
MT-PT4609236R	9	236
MT-PT4611276R	11	276
MT-PT4613316R	13	316



LPS Proximal Lateral Tibia Plate 4.5 | Material : TAN | Thickness : 4.5mm

PRODUCT CODE	HOLES	LENGTH
MT-PT4805140L	5	140
MT-PT4807180L	7	180
MT-PT4809220L	9	220
MT-PT4811260L	11	260
MT-PT4813300L	13	300
MT-PT4805140R	5	140
MT-PT4807180R	7	180
MT-PT4809220R	9	220
MT-PT4811260R	11	260
MT-PT4813300R	13	300



LPS T Plate 4.5mm | Material : TiCP | Thickness : 2.5mm

PRODUCT CODE	HOLES	LENGTH
MT-PT3003068	3	68
MT-PT3004084	4	84
MT-PT3005100	5	100
MT-PT3006116	6	116
MT-PT3007132	7	132
MT-PT3008148	8	148
MT-PT3010180	10	180
MT-PT3012212	12	212



ARMAR™ 4.5mm Meril Locking Plate System

LPS T Buttress Plate 4.5 | Material : TiCP | Thickness : 2.5mm

PRODUCT CODE	HOLES	LENGTH
MT-PT2804083	4	83
MT-PT2805099	5	99
MT-PT2806115	6	115



LPS L Buttress Plate 4.5 | Material : TiCP | Thickness : 2.5mm

PRODUCT CODE	HOLES	LENGTH
MT-PT2703069L	3	72
MT-PT2704085L	4	88
MT-PT2705101L	5	104
MT-PT2706117L	6	120
MT-PT2703069R	3	72
MT-PT2704085R	4	88
MT-PT2705101R	5	104
MT-PT2706117R	6	120



LPS Large Metaphyseal Plate 4.5 | Material : TiCP | Thickness : 4.2mm


PRODUCT CODE	HOLES	LENGTH
MT-PT0403118	8	118
MT-PT0404136	9	136
MT-PT0405154	10	154
MT-PT0406172	11	172
MT-PT0407190	12	190
MT-PT0408208	13	208
MT-PT0409226	14	226
MT-PT0411262	16	262
MT-PT0413298	18	298
MT-PT0415334	20	334



MBOSS™ 4.5mm Screw

Ø 5.0mm Locking Screw, Self Tapping | Material : TiCP
Head Diameter : 7.0mm | Outer Diameter : 5.0mm | Core Diameter : 4.4mm | Pitch : 1.0mm

PRODUCT CODE	LENGTH
MT-ST2850014	14
MT-ST2850016	16
MT-ST2850018	18
MT-ST2850020	20
MT-ST2850022	22
MT-ST2850024	24
MT-ST2850026	26
MT-ST2850028	28
MT-ST2850030	30
MT-ST2850032	32
MT-ST2850034	34
MT-ST2850036	36
MT-ST2850038	38
MT-ST2850040	40
MT-ST2850042	42
MT-ST2850044	44
MT-ST2850046	46
MT-ST2850048	48
MT-ST2850050	50
MT-ST2850055	55
MT-ST2850060	60
MT-ST2850065	65
MT-ST2850070	70
MT-ST2850075	75
MT-ST2850080	80
MT-ST2850085	85
MT-ST2850090	90



MBOSS™ 4.5mm Screw

Ø4.5mm Cortex Screw, Self Tapping | Material : TiCP
Head Diameter : 8.0mm | Outer Diameter : 4.5mm | Core Diameter : 3.0mm | Pitch : 1.75mm

PRODUCT CODE

LENGTH

MT-ST2445014	14
MT-ST2445016	16
MT-ST2445018	18
MT-ST2445020	20
MT-ST2445022	22
MT-ST2445024	24
MT-ST2445026	26
MT-ST2445028	28
MT-ST2445030	30
MT-ST2445032	32
MT-ST2445034	34
MT-ST2445036	36
MT-ST2445038	38
MT-ST2445040	40
MT-ST2445042	42
MT-ST2445044	44
MT-ST2445046	46
MT-ST2445048	48
MT-ST2445050	50
MT-ST2445052	52
MT-ST2445054	54
MT-ST2445056	56
MT-ST2445058	58
MT-ST2445060	60
MT-ST2445062	62
MT-ST2445064	64
MT-ST2445066	66
MT-ST2445068	68
MT-ST2445070	70
MT-ST2445072	72
MT-ST2445074	74
MT-ST2445076	76
MT-ST2445078	78
MT-ST2445080	80
MT-ST2445085	85
MT-ST2445090	90
MT-ST2445095	95
MT-ST2445100	100



MBOSS™ 4.5mm Screw

Ø6.5mm Cancellous Screw, Full Thread Length | Material : TiCP
Head Diameter : 8.0mm | Outer Diameter : 6.5mm | Core Diameter : 3.2mm | Pitch : 2.75mm

PRODUCT CODE

LENGTH

MT-ST3665020	20
MT-ST3665025	25
MT-ST3665030	30
MT-ST3665035	35
MT-ST3665040	40
MT-ST3665045	45
MT-ST3665050	50
MT-ST3665055	55
MT-ST3665060	60
MT-ST3665065	65
MT-ST3665070	70
MT-ST3665075	75
MT-ST3665080	80
MT-ST3665085	85
MT-ST3665090	90
MT-ST3665095	95
MT-ST3665100	100
MT-ST3665105	105
MT-ST3665110	110



MBOSS™ 4.5mm Screw

Ø 6.5mm Cancellous Screw, 16mm Thread Length | Material : TiCP
Head Diameter : 8.0mm | Outer Diameter : 6.5mm | Core Diameter : 3.2mm | Pitch : 2.75mm

PRODUCT CODE	LENGTH
MT-ST3165030	30
MT-ST3165035	35
MT-ST3165040	40
MT-ST3165045	45
MT-ST3165050	50
MT-ST3165055	55
MT-ST3165060	60
MT-ST3165065	65
MT-ST3165070	70
MT-ST3165075	75
MT-ST3165080	80
MT-ST3165085	85
MT-ST3165090	90
MT-ST3165095	95
MT-ST3165100	100
MT-ST3165105	105
MT-ST3165110	110
MT-ST3165115	115
MT-ST3165120	120



MBOSS™ 4.5mm Screw

Ø 6.5mm Cancellous Screw, 32mm Thread Length | Material : TiCP
Head Diameter : 8.0mm | Outer Diameter : 6.5mm | Core Diameter : 3.2mm | Pitch : 2.75mm

PRODUCT CODE	LENGTH
MT-ST3265045	45
MT-ST3265050	50
MT-ST3265055	55
MT-ST3265060	60
MT-ST3265065	65
MT-ST3265070	70
MT-ST3265075	75
MT-ST3265080	80
MT-ST3265085	85
MT-ST3265090	90
MT-ST3265095	95
MT-ST3265100	100
MT-ST3265105	105
MT-ST3265110	110
MT-ST3265115	115
MT-ST3265120	120
MT-ST3265125	125
MT-ST3265130	130
MT-ST3265135	135
MT-ST3265140	140
MT-ST3265145	145
MT-ST3265150	150




MBOSS™ 4.5mm Screw

Ø 5.0mm Locking Cancellous Screw, Self Tapping | Material : TiCP
Head Diameter : 7.0 mm | Outer Diameter : 5.0 mm | Core Diameter : 3.2 mm | Pitch : 2.75 mm

PRODUCT CODE	LENGTH
MT-ST8550050	50
MT-ST8550055	55
MT-ST8550060	60
MT-ST8550065	65
MT-ST8550070	70
MT-ST8550075	75
MT-ST8550080	80
MT-ST8550085	85
MT-ST8550090	90



ARMAR™ 4.5mm Locking Plate System Instrument Set



ARMAR™ 4.5mm Locking Plate System Instrument Set

Part No:MT-INLLPS001	Description
	Drill bit Ø3.2mm, Length 145mm, 2-flute, for Quick Coupling
	
Part No:MT-INLLPS002	Description
	Drill bit Ø4.3mm, Length 270 mm, 2-flute, for Quick Coupling
	
Part No:MT-INLLPS003	Description
	Drill bit Ø4.5mm, Length 145mm, 2-flute, for Quick Coupling
	
Part No:MT-INLLPS005	Description
	Tap for Cortex Screws Ø4.5mm, Length 130mm for quick coupling
	
Part No:MT-INLLPS015	Description
	Ø3.2 Insert Drill Sleeve
	
Part No:MT-INLLPS010	Description
	LPS Drill Sleeve 5.0, for Drill Bits Ø4.3mm
	

ARMAR™ 4.5mm Locking Plate System Instrument Set

Part No:MT-INLLPS017

Description

Drill Sleeve OD4.3mm & ID2.0mm



Part No:MT-INLLPS013

Description

Double Drill Guide 4.5mm & 3.2mm, for Neutral and Load position



Part No:MT-INLLPS011

Description

Universal Drill Guide 4.5mm & 3.2mm



Part No:MT-INLLPS024

Description

Screwdriver, Hexagonal, A/F 3.5mm, with Groove & Holding Sleeve



Part No:MT-INLLPS021

Description

Screwdriver, hexagonal, large, Ø3.5mm, with Groove, Length 305mm



ARMAR™ 4.5mm Locking Plate System Instrument Set

Part No:MT-INLLPS019

Description

Torque Limiter Screw Driver 3.5mm Hex Self Holding 4Nm



Part No:MT-INSLPS026

Description

Wrench-SLPS



Part No:MT-INSLPS031

Description

T-Handle with Quick coupling SLPS



Part No:MT-INSLPS017

Description

Guide Wire - Krischner - SLPS



Part No:MT-INLLPS025

Description

Depth Gauge for Screws 4.5/6.5mm, measuring range up to 100mm



ARMAR™ 4.5mm Locking Plate System Instrument Set

Part No:MT-INLLPS012

Description


Double Drill Guide 4.5mm & 3.2mm



Part No:MT-INLLPS039

Description

Instrument Case with Lid - LLPS



MBOSS™ 4.5mm Cannulated Screw

4.5mm Cannulated Screw, Self drilling, Short Thread

Material : TAN | Head Diameter : 6.5 | Outer Diameter : 4.5 | Core Diameter : 2.7 | Pitch : 1.75

PRODUCT CODE	LENGTH
--------------	--------

MT-ST2345020	20
MT-ST2345022	22
MT-ST2345024	24
MT-ST2345026	26
MT-ST2345028	28
MT-ST2345030	30
MT-ST2345032	32
MT-ST2345034	34
MT-ST2345036	36
MT-ST2345038	38
MT-ST2345040	40
MT-ST2345042	42
MT-ST2345044	44
MT-ST2345046	46
MT-ST2345048	48
MT-ST2345050	50
MT-ST2345052	52
MT-ST2345054	54
MT-ST2345056	56
MT-ST2345060	60
MT-ST2345064	64
MT-ST2345068	68
MT-ST2345072	72
MT-ST2345076	76
MT-ST2345080	80



MBOSS™ 4.5mm Cannulated Screw

Ø10.0mm Washer

Material : TAV | Diameter : 10.0mm | Thickness : 1.2

PRODUCT CODE	LENGTH
--------------	--------

MT-WT4210000

-



MBOSS™ 4.5 Cancellous Cannulated Screw Instrument Set



MBOSS™ 4.5mm Cancellous Cannulated Screw Instrument set

Part No:MT-INCS45001

Description

Drill bit Ø3.2mm, with 1.75mm Cannulation Length 170/140mm



Part No:MT-INCS45002

Description

Cannulated Countersink for Cannulated Ø4.5mm Screws



Part No:MT-INCS45003

Description

Parallel Guide for Guide wire Ø1.6mm, Adjustable for Cannulated Screw



Part No:MT-INCS45005

Description

Drill sleeve 3.2/1.6, for percutaneous insertion of cannulated Ø4.5mm Screws



MBOSS™ 4.5mm Cancellous Cannulated Screw Instrument Set

Part No:MT-INCS45006

Description

Drill sleeve 7/3.2, for percutaneous insertion of cannulated Ø4.5mm Screws



Part No:MT-INCS45007

Description

Protection sleeve 9.5/7, percutaneous insertion of cannulated Ø4.5mm Screws



Part No:MT-INCS45008

Description

Hexagonal Cannulated Screwdriver for cannulated Ø4.5mm Screws



Part No:MT-INCS45009

Description

Hexagonal Screwdriver for Ø4.5mm Screws



Part No:MT-INCS45010

Description

Direct measuring device for cannulated Ø4.5mm Screws



MBOSS™ 4.5mm Cancellous Cannulated Screw Instrument Set

Part No:MT-INCS45013

Description

Cleaning Pin



Part No:MT-INCS45019

Description

Guide Wire φ1.5×150 - 1



Part No:MT-INCS45020

Description

Guide Wire φ1.5×150 - 2



Part No:MT-INCS45021

Description


Trocar φ 1.5mm, for percutaneous insertion of cannulated Ø4.5mm Screws



Part No:MT-INCS45018

Description

Instrument Case With Lid



MBOSS™ 6.5mm Cannulated Screw

6.5mm Cannulated Screw, Self drilling, 16mm Thread Length

Material : TAN | Head Diameter : 8.0 | Outer Diameter : 6.5 | Core Diameter : 4.5 | Pitch : 2.75

PRODUCT CODE	LENGTH
--------------	--------

MT-ST3965030	30
MT-ST3965035	35
MT-ST3965040	40
MT-ST3965045	45
MT-ST3965050	50
MT-ST3965055	55
MT-ST3965060	60
MT-ST3965065	65
MT-ST3965070	70
MT-ST3965075	75
MT-ST3965080	80
MT-ST3965085	85
MT-ST3965090	90
MT-ST3965095	95
MT-ST3965100	100
MT-ST3965105	105
MT-ST3965110	110



6.5mm Cannulated Screw, Self drilling, 32mm Thread Length

Material : TAN | Head Diameter : 8.0 | Outer Diameter : 6.5 | Core Diameter : 4.5 | Pitch : 2.75

PRODUCT CODE	LENGTH
--------------	--------

MT-ST3765045	45
MT-ST3765050	50
MT-ST3765055	55
MT-ST3765060	60
MT-ST3765065	65
MT-ST3765070	70
MT-ST3765075	75
MT-ST3765080	80
MT-ST3765085	85
MT-ST3765090	90
MT-ST3765095	95
MT-ST3765100	100
MT-ST3965105	105
MT-ST3965110	110



MBOSS™ 6.5mm Cannulated Screw

Ø13.0mm Washer

Material : TAV | Diameter : 13.0mm | Thickness : 1.5

PRODUCT CODE	LENGTH
--------------	--------

MT-WT4213000

-



MBOSS™ 6.5mm Cancellous Cannulated Screw Instrument Set



MBOSS™ 6.5mm Cancellous Cannulated Screw Instrument Set

Part No:MT-INCC65001

Description


5.0mm Cannulated Drill Bit 300mm



Part No:MT-INCC65002

Description

Screw Measuring Device for guide wire



Part No:MT-INCC65003

Description

2.5mm x 300mm Threaded Guide Wire Trocar point



Part No:MT-INCC65005

Description

Screwdriver hexagonal cannulated 4.0mm



Part No:MT-INCC65006

Description


Fixation Sleeve for Drill Bits-CC65



Part No:MT-INCC65011

Description

Cleaning brush for Cc65



MBOSS™ 6.5mm Cancellous Cannulated Screw Instrument Set

Part No:MT-INCC65012

Description

12.0mm / 8.5mm Protection Sleeve



Part No:MT-INCC65013

Description


Drill Sleeve 8.5mm / 2.8mm for Guide Wire Ø2.5mm



Part No:MT-INCC65015

Description

Cleaning Pin for Cc65



Part No:MT-INCC65016

Description

15.5mm / 13.0mm Washer Sleeve



Part No:MT-INCC65017

Description

Double Drill Sleeve for guide wire



MBOSS™ 6.5mm Cancellous Cannulated Screw Instrument Set

Part No:MT-INSLPS033

Description

Screw Holding Forceps



Part No:MT-INCC65007

Description

Wrench for Fixation Sleeve-CC65



Part No:MT-INCC65009

Description

Screwdriver, hexagonal 4.0mm



Part No:MT-INCC65020

Description

2.5mm Trocar



Part No:MT-INCC65021

Description

Drill Sleeve 8.5mm / 5.5mm for Ø5.0mm Dril Bit



MBOSS™ 6.5mm Cancellous Cannulated Screw Instrument Set

Part No:MT-INCC65014

Description

Instrument Case with Lid - CS65



ARMAR™ Proximal Humeral Locking Plate System Instrument Set



ARMAR™ Proximal Humeral Locking Plate System Instrument Set

Part No:MT-INPHLP001

Description

Aiming Device for Proximal Humeral Plate



Part No:MT-INPHLP002

Description

3.5 LPS Drill Sleeve for Ø2.8mm Drill Bits



Part No:MT-INPHLP003

Description

Centering Sleeve for Aiming Device



Part No:MT-INPHLP004

Description

Drill Sleeve for Aiming Device



Part No:MT-INPHLP005

Description

Drill Sleeve for Guide Wire Ø1.5mm



ARMAR™ Proximal Humeral Locking Plate System Instrument Set

Part No:MT-INPHLP006

Description

Depth Gauge for PHLP



Part No:MT-INPHLP007

Description

Guide Wire Ø1.5, length 150mm



Part No:MT-INPHLP008

Description

Instrument Case with Lid - PHLP



ARMAR™ Proximal Lateral Tibia Locking Plate System Instrument Set



ARMAR™ Proximal Lateral Tibia Locking Plate System Instrument Set

Part No:MT-INPLTP001

Description

Drill Sleeve for Insertion guide-PLTP



Part No:MT-INPLTP002

Description

Drill Sleeve for Ø4.3mm Drill bit-PLTP



Part No:MT-INPLTP003

Description


Fixation Bolt for Insetion Guide



Part No:MT-INPLTP004

Description

Trocar for PLTP



Part No:MT-INPLTP005

Description

Depth Gauge for PLTP



Part No:MT-INPLTP006

Description

Centering Sleeve for Guide Wire Ø2.5mm



ARMAR™ Proximal Lateral Tibia Locking Plate System Instrument Set

Part No:MT-INPLTP007

Description

Insertion Guide for Proximal Tibia Right



Part No:MT-INPLTP008

Description

Insertion Guide for Proximal Tibia Left



Part No:MT-INDFLP009

Description


Guide Wire Ø2.5, length 250mm



Part No:MT-INDFLP010

Description

Wrench for Fixation Sleeve



ARMAR™ Proximal Lateral Tibia Locking Plate System Instrument Set

Part No:MT-INDFLP011

Description

Fixation Sleeve for Drill Bits



Part No:MT-INDFLP012

Description

Pulling Device



Part No:MT-INDFLP013

Description

Wrench for Fixation Bolt



Part No:MT-INDFLP014

Description

Stopper for Insertion Guide



Part No:MT-INPLTP015

Description

Screwdriver Shaft A/F3.5, hexagonal Large



ARMAR™ Proximal Lateral Tibia Locking Plate System Instrument Set

Part No:MT-INPLTP016

Description

Screwdriver, Hexagonal, Large, A/F3.5mm



Part No:MT-INPLTP018

Description

Drill Bit Ø4.3mm, Length 300mm, 2 flute, for Quick coupling



Part No:MT-INPLTP019

Description

Instrument Case with Lid - PLTP



ARMAR™ Distal Femoral Locking Plate System Instrument Set



ARMAR™ Distal Femoral Locking Plate System Instrument Set

Part No:MT-INDFLP001

Description

Insertion Guide for Distal Femur, Left



Part No:MT-INDFLP002

Description

Centering Sleeve for Guide Wire Ø2.5mm-DFLP



Part No:MT-INDFLP003

Description

Drill Sleeve for Insertion guide-DFLP



Part No:MT-INDFLP004

Description


Fixation Bolt for Insertion Guide



Part No:MT-INDFLP005

Description

Drill Sleeve for Ø4.3mm Drill bit-DFLP



ARMAR™ Distal Femoral Locking Plate System Instrument Set

Part No:MT-INDFLP006

Description

Trocar for DFLP
 

Part No:MT-INDFLP007

Description

Depth Gauge for DFLP
 

Part No:MT-INDFLP008

Description

Insertion Guide for Distal Femur, Right
 

Part No:MT-INDFLP009

Description

Guide Wire Ø2.5, length 250mm
 

Part No:MT-INDFLP010

Description

Wrench for Fixation Sleeve
 

ARMAR™ Distal Femoral Locking Plate System Instrument Set

Part No:MT-INDFLP011

Description

Fixation Sleeve for Drill Bits
 

Part No:MT-INDFLP012

Description

Pulling Device
 

Part No:MT-INDFLP013

Description

Wrench for Fixation Bolt
 

Part No:MT-INDFLP014

Description

Stopper for Insertion Guide
 

Part No:MT-INDFLP015

Description

Screwdriver Shaft A/F3.5, hexagonal Large, self-holding sleeve
 

ARMAR™ Distal Femoral Locking Plate System Instrument Set

Part No:MT-INDFLP016

Description

Screwdriver, Hexagonal, Large, A/F3.5mm, with Groove



Part No:MT-INDFLP018

Description

Drill Bit Ø4.3mm, Length 300mm, 2 flute, for Quick coupling



Part No:MT-INDFLP019

Description

Instrument Case with Lid - DFLP

