

Corporate Headquarters, Vapi

Meril Healthcare Pvt. Ltd.
H1-H3, Meril Park
Survey No. 135/2/B & 174/2
Muktanand Marg, Chala
Vapi 396 191. Gujarat. India.
T: +91 260 3052 100
F: +91 260 3052 125
E: orthopedics@merillife.com

International Sales & Marketing

Meril Healthcare Pvt. Ltd.
612, B-Wing, Bonanaza
Sahar Plaza, Andheri East
Mumbai 400 059
Maharashtra. India.
T: +91 22 3935 0700
F: +91 22 4047 9717

India Sales & Marketing

Meril Healthcare Pvt. Ltd.
512, Midas, Sahar Plaza, J.B. Nagar
Andheri East, Mumbai 400 059
Maharashtra. India.
T: +91 22 4047 9797

Meril South America

DOC MED LTDA
1079 – Cep: 04077-003 - Moema
Sao Paulo, Brazil
T/F: +55 11 3624 5935/6

Meril GmbH

Bornheimer Strasse 135 – 137
D – 53119 Bonn, Germany
T/F: +49 228 7100 4000/1

Meril Tibbi Cihazlar

Imalat Ve Ticaret A.S.
6, Mimar Sinan Mah.,
Cavusbasi Cad, Ozde Sok
Aydin Eksi Is Merkezi Kat:1
Cekmekoy/Istanbul, Turkey
T: +90 53 2272 5172

Meril

Healthcare



KET™ Catalogue

SS Plate | Screw | Instrument

Cardiovascular

Orthopedics

Diagnostics

Endo-Surgery

ENT

Meril

More to Life

More to Life

Meril is a global medical device company dedicated to the design and development of novel, clinically relevant, state-of-the-art and best-in-class devices to alleviate human suffering and improve the quality of life. We span a broad operational canvas ranging from vascular interventional devices to orthopedics, in-vitro diagnostics and endo-surgery.

We share an enduring perseverance to advance healthcare solutions, for more patients to live longer and healthier lives. We have a strong commitment towards R&D and adhere to the best quality standards in manufacturing, scientific communication and distribution known today.

Established in 2006, Meril's launch was in line with the healthcare diversification plan of Bilakhia Group. Located in Vapi, 150km north of Mumbai, and ensconced within the serene hill-sides of Chala town, Meril sits proud in its ultra modern manufacturing facilities.

Where research becomes the basis for discovering new ways of alleviating human suffering and empowering an individual with good health

ri Research
Empowering
Individual




Infrastructure

Situated on 180 acres of land and built over a 300,000 sq.ft. area, Meril's ISO 13485 & cGMP certified ultra-modern manufacturing facility is constructed to satisfy the stringent needs of our cardiovascular, orthopedics, diagnostics and endo-surgery businesses, as well as our R&D endeavors.

Fully integrated manufacturing systems ensure backward integration, flawless man-material movement and complete control over processes to eliminate production errors. All manufacturing and sterilization processes are conducted in-house, in addition to the analytical and microbiological QA/QC tests required to meet the world-class production standards.



Meril's 3,200+ personnel are a strong, experienced team comprising designers, engineers, chemists, microbiologists, regulatory affair experts, R&D scientists, clinical affairs experts, legal, finance, sales and marketing professionals with innovative capabilities. The team is continuously striving towards improving and saving lives every day by providing revolutionary solutions for diagnosis, prevention and treatment.



Global Market Presence

Currently Meril has direct sales & representatives offices in India, Germany, Turkey, USA, China and Brazil and strong distribution model spread over 102 countries across continents.



Meril Orthopedic

Meril is known for an uncompromising attitude when it comes to quality of design. Providing the medical world with zero defect devices is our only motto. Our orthopaedic trauma products Meril Locking Plate System 3.5, Meril Locking Plate System 4.5, Meril Compression Plate System 3.5, Meril Compression Plate System 4.5, Meril Dynamic Hip Fixation System, and Meril Locking Distal Radial System 2.4 are specifically designed to treat fractures, bone tumour and deformities. They boast of high fatigue strength, proper fixation and negligible soft tissue damage.

We have been providing trauma implants for the foot, shoulder, hip, hand, pelvis, and condyles. In order to ensure smooth and trouble-free operating procedures leading medical professionals rely only on our Trauma products. We have carried out performance testing to prove our design for its endurance and durability and the outcome of these tests prove that they are going to create a dynamic image in the history of joint replacement.



Index

• KET™ 2.4mm Meril Locking Plate System	1
• KET™ 3.5mm Meril Locking Plate System	13
• KET™ 3.5mm Small Fragment Plate System	18
• KET™ Elbow-Clavicle Plate System	31
• KET™ 4.5mm Locking Plate System	45
• KET™ 4.5mm Large Fragment Plate System	48
• KET™ DHF & DSP Plates System	57

KET™ 2.4mm Meril Locking Plate System

LPS Distal Radial Plate 2.4, Volar
Material : SS 316 LVM | Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PS5303043L	5	3	43
MT-PS5305061L	5	5	61
MT-PS5303043R	5	3	43
MT-PS5305061R	5	5	61



LPS Distal Radial Plate 2.4, Volar Buttress
Material : SS 316 LVM | Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PS5304057L	5	3	52
MT-PS5304057R	5	3	52



LPS Distal Radial Plate 2.4, Extra-articular
Material : SS 316 LVM | Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PS5303048L	5	3	48
MT-PS5305066L	5	5	66
MT-PS5304047L	4	3	47
MT-PS5304065L	4	5	65
MT-PS5303048R	5	3	48
MT-PS5305066R	5	5	66
MT-PS5304047R	4	3	47
MT-PS5304065R	4	5	65



LPS T-Distal Radial Plate 2.4
Material : SS 316 LVM | Thickness : 1.6mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PS4203040	3	3	40
MT-PS4203049	3	4	49



KET™ 2.4mm Meril Locking Plate System


LPS Distal Radial Plate 2.4, Straight
Material : SS 316 LVM | Thickness : 1.6mm

PRODUCT CODE	HOLES	LENGTH
MT-PS4205048	5	48
MT-PS4206057	6	57



LPS L-Distal Radial Plate 2.4 Right Angled
Material : SS 316 LVM | Thickness : 1.6mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PS4202040L	2	3	40
MT-PS4202049L	2	4	49
MT-PS4203040L	3	3	40
MT-PS4203049L	3	4	49
MT-PS4202040R	2	3	40
MT-PS4202049R	2	4	49
MT-PS4203040R	3	3	40
MT-PS4203049R	3	4	49



LPS L-Distal Radial Plate 2.4, Oblique Angled
Material : SS 316 LVM | Thickness : 1.6mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PS4203043L	3	3	43
MT-PS4204052L	3	4	52
MT-PS4203043R	3	3	43
MT-PS4204052R	3	4	52



LPS Distal Radial Plate 2.4,Long
Material : SS 316 LVM | Thickness : 1.6mm


PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PS5308100	4	8	100
MT-PS5310120	4	10	120
MT-PS5312140	4	12	140



KET™ 2.4mm Meril Locking Plate System

2.4mm LPS Volar Rim Distal Radial Plate | Material : SS 316 LVM | Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PS7106057R	6	5	57
MT-PS7106057L	6	5	57
MT-PS7107057R	7	5	57
MT-PS7107057L	7	5	57



2.4mm LPS Double-Column Distal Radial Plate, Volar 19.5mm | Material : SS 316 LVM
Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PS6802042R	6	2	42
MT-PS6803051R	6	3	51
MT-PS6804063R	6	4	63
MT-PS6805072R	6	5	72
MT-PS6802042L	6	2	42
MT-PS6803051L	6	3	51
MT-PS6804063L	6	4	63
MT-PS6805072L	6	5	72



2.4mm LPS Double-Column Distal Radial Plate, Volar 22.0mm | Material : SS 316 LVM
Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PS6902045R	6	2	45
MT-PS6903054R	6	3	54
MT-PS6904066R	6	4	66
MT-PS6905075R	6	5	75
MT-PS6902045L	6	2	45
MT-PS6903054L	6	3	54
MT-PS6904066L	6	4	66
MT-PS6905075L	6	5	75



KET™ 2.4mm Meril Locking Plate System

2.4mm LPS Double-Column Distal Radial Plate, Volar 25.5mm | Material : SS 316 LVM
Thickness : 2.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PS7002047R	7	2	47
MT-PS7003056R	7	3	56
MT-PS7004068R	7	4	68
MT-PS7005077R	7	5	77
MT-PS7002047L	7	2	47
MT-PS7003056L	7	3	56
MT-PS7004068L	7	4	68
MT-PS7005077L	7	5	77



KET™ 2.4mm Screw

2.4mm Cortex Screw, T8 Self Tapping | Material : SS 316 LVM | Head Diameter : 4.0mm
Outer Diameter : 2.4mm | Core Diameter : 1.7mm | Pitch : 1mm

PRODUCT CODE

LENGTH

MT-SS0224008	8
MT-SS0224010	10
MT-SS0224012	12
MT-SS0224014	14
MT-SS0224016	16
MT-SS0224018	18
MT-SS0224020	20
MT-SS0224022	22
MT-SS0224024	24
MT-SS0224026	26
MT-SS0224028	28
MT-SS0224030	30




2.7mm Cortex Screw, T8 Self Tapping

Material : SS 316 LVM | Head Diameter : 5.0mm
Outer Diameter : 2.7mm | Core Diameter : 1.9mm | Pitch : 1mm

PRODUCT CODE

LENGTH

MT-SS0527010	10
MT-SS0527012	12
MT-SS0527014	14
MT-SS0527016	16
MT-SS0527018	18
MT-SS0527020	20
MT-SS0527022	22
MT-SS0527024	24
MT-SS0527026	26
MT-SS0527028	28
MT-SS0527030	30



KET™ 2.4mm Screw

2.4mm Locking Screw, T8 Self Tapping | Material : SS 316 LVM | Head Diameter : 3.6mm
Outer Diameter : 2.4mm | Core Diameter : 1.9mm | Pitch : 0.6mm

PRODUCT CODE

LENGTH

MT-SS0124008	8
MT-SS0124010	10
MT-SS0124012	12
MT-SS0124014	14
MT-SS0124016	16
MT-SS0124018	18
MT-SS0124020	20
MT-SS0124022	22
MT-SS0124024	24
MT-SS0124026	26
MT-SS0124028	28
MT-SS0124030	30



KET™ 2.4mm Screw

2.4mm Headless Compression Screw Short Thread

Material : SS 316 LVM | Head Diameter : 3.1 | Outer Diameter : 2.4 | Core Diameter : 1.9 | Pitch : 1.0

PRODUCT CODE

LENGTH

MT-SS0324009	09
MT-SS0324010	10
MT-SS0324012	12
MT-SS0324014	14
MT-SS0324016	16
MT-SS0324018	18
MT-SS0324020	20
MT-SS0324022	22
MT-SS0324024	24
MT-SS0324026	26
MT-SS0324028	28
MT-SS0324030	30
MT-SS0324032	32
MT-SS0324034	34
MT-SS0324036	36
MT-SS0324038	38
MT-SS0324040	40



2.4mm Headless Compression Screw Long Thread

Material : SS 316 LVM | Head Diameter : 3.1 | Outer Diameter : 2.4 | Core Diameter : 1.9 | Pitch : 1.0


PRODUCT CODE

LENGTH

MT-SS0424016	16
MT-SS0424018	18
MT-SS0424020	20
MT-SS0424022	22
MT-SS0424024	24
MT-SS0424026	26
MT-SS0424028	28
MT-SS0424030	30



KET™ 2.4mm Locking Plate system Instrument Set



KET™ 2.4mm Locking Plate system Instrument Set

Part No:MT-INDR24001

Description

Drill bit Ø1.8mm, Length 110mm, 2-flute, for Quick Coupling



Part No:MT-INDR24002

Description

Drill Bit Ø2.0mm, Length 110mm, 2-flute, for Quick Coupling



Part No:MT-INDR24008

Description

Universal Drill Guide 2.7



Part No:MT-INDR24007

Description

Universal Drill Guide 2.4



Part No:MT-INDR24005

Description

LPS Drill Sleeve 2.4, with Scale up to 30mm, for Drill bits Ø1.8mm



KET™ 2.4mm Locking Plate system Instrument Set

Part No:MT-INDR24006

Description

LPS Drill Sleeve 2.7 with Scale up to 30mm, for Drill bits Ø2.0mm



Part No:MT-INDR24015

Description

Screwdriver shaft, T8, self-holding



Part No:MT-INDR24013

Description

Manual 0.8Nm Torque Limiter



Part No:MT-INDR24020

Description

Quick Change Handle Large with Quick Coupling, Length 138



Part No:MT-INDR24016

Description

Screw Driver Stardrive T8 Self Holding



Part No:MT-INDR24018

Description

Depth Gauge for Screws 2.0/2.4mm, measuring range up to 40mm



KET™ 2.4mm Locking Plate system Instrument Set

Part No:MT-INDR24021

Description

Grasping sets of 2.4/2.7 locking screw, Screw Driver T8



Part No:MT-INDR24028

Description


Ø1.25mm Guide Wire Length 150mm



Part No:MT-INDR24029

Description

Instrument Case with Lid



KET™ 3.5mm Meril Locking Plate System

3.5mm LPS-Small Compression Plate
Material : SS 316 LVM | Thickness : 3.3mm

PRODUCT CODE	HOLES	LENGTH
MT-PS0804059	4	59
MT-PS0805072	5	72
MT-PS0806085	6	85
MT-PS0807098	7	98
MT-PS0808111	8	111
MT-PS0809124	9	124
MT-PS0810137	10	137
MT-PS0811150	11	150
MT-PS0812163	12	163



3.5mm LPS-Small Metaphyseal Plate
Material : SS 316 LVM | Thickness : 3.3mm

PRODUCT CODE	HOLES	LENGTH
MT-PS1606086	6	86
MT-PS1607099	7	99
MT-PS1608112	8	112
MT-PS1609125	9	125
MT-PS1610138	10	138



3.5mm LPS-Reconstruction Round Hole Plate
Material : SS 316 LVM | Thickness : 3.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS0704046	4	46
MT-PS0705058	5	58
MT-PS0706070	6	70
MT-PS0707082	7	82
MT-PS0708094	8	94
MT-PS0709106	9	106
MT-PS0710118	10	118
MT-PS0712142	12	142



KET™ 3.5mm Meril Locking Plate System

3.5mm LPS-Reconstruction Combined Hole Plate
Material : SS 316 LVM | Thickness : 3.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS1804056	4	56
MT-PS1805070	5	70
MT-PS1806084	6	84
MT-PS1807098	7	98
MT-PS1808112	8	112
MT-PS1809126	9	126
MT-PS1810140	10	140
MT-PS1812168	12	168



3.5mm LPS T Plate Right Angled
Material : SS 316 LVM | Thickness : 3.0mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PS1903050R	3	3	46
MT-PS1905067R	3	5	67
MT-PS1904056R	4	4	57
MT-PS1906078R	4	6	78



3.5mm LPS T Plate Oblique
Material : SS 316 LVM | Thickness : 1.6mm

PRODUCT CODE	HOLES	SHAFT HOLES	LENGTH
MT-PS2003052R	3	3	52
MT-PS2003052L	3	3	52
MT-PS2004063R	3	4	63
MT-PS2004063L	3	4	63
MT-PS2005074R	3	5	74
MT-PS2005074L	3	5	74



3.5mm LPS Cloverleaf Plate
Material : SS 316 LVM | Thickness : 2.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PT1303093	3	92.5
MT-PT1304107	4	107
MT-PT1305122	5	121.5
MT-PT1306136	6	136



KET™ 3.5mm Meril Locking Plate System

3.5mm LPS Proximal Humerus Plate, Short

Material : SS 316 LVM | Thickness : 3.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS5103090	3	90
MT-PS5105114	5	114



3.5mm LPS Proximal Humerus Plate, Long

Material : SS 316 LVM | Thickness : 4.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS5205142	5	142
MT-PS5206160	6	160
MT-PS5208196	8	196
MT-PS5210232	10	232



3.5mm LPS Proximal Tibia Plate

Material : SS 316 LVM | Thickness : 4.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS1704081R	4	76
MT-PS1706107R	6	102
MT-PS1708133R	8	128
MT-PS1710159R	10	154
MT-PS1712185R	12	180
MT-PS1704081L	4	76
MT-PS1706107L	6	102
MT-PS1708133L	8	128
MT-PS1710159L	10	154
MT-PS1712185L	12	180



KET™ 3.5mm Meril Locking Plate System

3.5mm LPS Medial Proximal Tibial Plate

Material : SS 316 LVM | Thickness : 4.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS1404093L	4	93
MT-PS1406119L	6	119
MT-PS1408145L	8	145
MT-PS1410171L	10	171
MT-PS1412197L	12	197
MT-PS1404093R	4	93
MT-PS1406119R	6	119
MT-PS1408145R	8	145
MT-PS1410171R	10	171
MT-PS1412197R	12	197



3.5mm LPS Medial Distal Tibia Low Bend Plate

Material : SS 316 LVM | Thickness : 4.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS4704109L	4	109
MT-PS4706135L	6	135
MT-PS4708161L	8	161
MT-PS4710187L	10	187
MT-PS4712213L	12	213
MT-PS4714239L	14	239
MT-PS4704109R	4	109
MT-PS4706135R	6	135
MT-PS4708161R	8	161
MT-PS4710187R	10	187
MT-PS4712213R	12	213
MT-PS4714239R	14	239



KET™ 3.5mm Meril Locking Plate System

3.5mm LPS One Third Tubular Plate

Material : SS 316 LVM | Thickness : 1.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS5005064	5	64
MT-PS5006076	6	76
MT-PS5007088	7	88
MT-PS5008100	8	100
MT-PS5009112	9	112



3.5mm LPS Anterior-Lateral Distal Tibial Plate

Material : SS 316 LVM | Thickness : 4.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS1205080L	5	80
MT-PS1207106L	7	106
MT-PS1209132L	9	132
MT-PS1211158L	11	158
MT-PS1213184L	13	184
MT-PS1205080R	5	80
MT-PS1207106R	7	106
MT-PS1209132R	9	132
MT-PS1211158R	11	158
MT-PS1213184R	13	184



3.5mm LPS Proximal Tibia Plate, Posteromedial

Material : SS 316 LVM | Thickness : 3.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PS7501069	1	69
MT-PS7502079	2	79
MT-PS7504105	4	105
MT-PS7506131	6	131
MT-PS7508157	8	157
MT-PS7510183	10	183



KET™ 3.5mm Meril Locking Plate System

3.5mm LPS Periarticular Proximal Lateral Tibia Plate

Material : SS 316 LVM | Thickness : 3.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PS7306112L	6	112
MT-PS7308138L	8	138
MT-PS7310164L	10	164
MT-PS7312190L	12	190
MT-PS7314216L	14	216
MT-PS7316242L	16	242
MT-PS7306112R	6	112
MT-PS7308138R	8	138
MT-PS7310164R	10	164
MT-PS7312190R	12	190
MT-PS7314216R	14	216
MT-PS7316242R	16	242



3.5mm LPS Periarticular Proximal Lateral Humeral Plate

Material : SS 316 LVM | Thickness : 3.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PS7404078L	4	78
MT-PS7406104L	6	104
MT-PS7408130L	8	130
MT-PS7412182L	12	182
MT-PS7416234L	16	234
MT-PS7404078R	4	78
MT-PS7406104R	6	104
MT-PS7408130R	8	130
MT-PS7412182R	12	182
MT-PS7416234R	16	234



KET™ 3.5mm Meril Locking Plate System

2.7mm/3.5mm LPS Lateral Distal Fibula Plate

Material : TiCP | Thickness : 3.0 mm

PRODUCT CODE	HOLES	LENGTH
MT-PS6704080R	4	80
MT-PS6706106R	6	106
MT-PS6708132R	8	132
MT-PS6710158R	10	158
MT-PS6704080L	4	80
MT-PS6706106L	6	106
MT-PS6708132L	8	132
MT-PS6710158L	10	158



KET™ 3.5mm Small Fragment Plate System

Small Compression Plate 3.5 | Material : SS 316 LVM | Thickness : 3.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS0602026	2	26
MT-PS0603038	3	38
MT-PS0604050	4	50
MT-PS0605062	5	62
MT-PS0606074	6	74
MT-PS0607086	7	86
MT-PS0608098	8	98
MT-PS0609110	9	110
MT-PS0610122	10	122
MT-PS0611134	11	134
MT-PS0612146	12	146



Low Contourable Dual Compression Plate 3.5 | Material : SS 316 LVM | Thickness : 3.3mm

PRODUCT CODE	HOLES	LENGTH
MT-PS2504051	4	51
MT-PS2505064	5	64
MT-PS2506077	6	77
MT-PS2507090	7	90
MT-PS2508103	8	103
MT-PS2509116	9	116
MT-PS2510129	10	129
MT-PS2511142	11	142
MT-PS2512155	12	155



KET™ 3.5mm Screw

3.5mm Cortex Screw, Self Tapping | Material : SS 316 LVM | Head Diameter : 6.0mm
Outer Diameter : 3.5mm | Inner Diameter : 2.4mm | Pitch : 1.75mm

PRODUCT CODE

LENGTH

MT-SS1035010	10
MT-SS1035012	12
MT-SS1035014	14
MT-SS1035016	16
MT-SS1035018	18
MT-SS1035020	20
MT-SS1035022	22
MT-SS1035024	24
MT-SS1035026	26
MT-SS1035028	28
MT-SS1035030	30
MT-SS1035032	32
MT-SS1035034	34
MT-SS1035036	36
MT-SS1035038	38
MT-SS1035040	40
MT-SS1035042	42
MT-SS1035044	44
MT-SS1035045	45
MT-SS1035046	46
MT-SS1035048	48
MT-SS1035050	50
MT-SS1035055	55
MT-SS1035060	60
MT-SS1035065	65
MT-SS1035070	70
MT-SS1035075	75
MT-SS1035080	80
MT-SS1035085	85
MT-SS1035090	90
MT-SS1035095	95
MT-SS1035100	100
MT-SS1035105	105
MT-SS1035110	110




KET™ 3.5mm Screw

3.5mm Locking Screw, Self Tapping | Material : SS 316 LVM | Head Diameter : 5.0mm
Outer Diameter : 3.5mm | Inner Diameter : 2.9mm | Pitch : 0.8mm

PRODUCT CODE

LENGTH

MT-SS1235010	10
MT-SS1235012	12
MT-SS1235014	14
MT-SS1235016	16
MT-SS1235018	18
MT-SS1235020	20
MT-SS1235022	22
MT-SS1235024	24
MT-SS1235026	26
MT-SS1235028	28
MT-SS1235030	30
MT-SS1235032	32
MT-SS1235034	34
MT-SS1235036	35
MT-SS1235038	36
MT-SS1235040	38
MT-SS1235042	40
MT-SS1235044	42
MT-SS1235046	46
MT-SS1235048	48
MT-SS1235050	50
MT-SS1235052	52
MT-SS1235054	54
MT-SS1235056	56
MT-SS1235058	58
MT-SS1235060	60
MT-SS1235065	65
MT-SS1235070	70
MT-SS1235075	75
MT-SS1235080	80
MT-SS1235085	85



KET™ 3.5mm Screw

3.5mm Locking Cancellous Screw, Self Tapping | Material : SS 316 LVM
Head Diameter : 5.0mm | Outer Diameter : 3.5mm | Inner Diameter : 2.0mm | Pitch : 0.8mm

PRODUCT CODE

LENGTH

MT-SS8435030	30
MT-SS8435050	50
MT-SS8435065	65
MT-SS8435070	70
MT-SS8435075	75
MT-SS8435080	80
MT-SS8435085	85
MT-SS8435090	90



4.0 mm Cancellous Screw, Full Thread | Material : SS 316 LVM | Head Diameter : 6.0mm
Outer Diameter : 4.0mm | Inner Diameter : 1.9mm | Pitch : 1.75mm

PRODUCT CODE

LENGTH

MT-SS1640010	10
MT-SS1640012	12
MT-SS1640014	14
MT-SS1640016	16
MT-SS1640018	18
MT-SS1640020	20
MT-SS1640022	22
MT-SS1640024	24
MT-SS1640026	26
MT-SS1640028	28
MT-SS1640030	30
MT-SS1640035	35
MT-SS1640040	40
MT-SS1640045	45
MT-SS1640050	50
MT-SS1640055	55
MT-SS1640060	60
MT-SS1640065	65
MT-SS1640070	70
MT-SS1640075	75
MT-SS1640080	80



KET™ 3.5mm Screw

4.0mm Cancellous Screw, Short Thread | Material : SS 316 LVM | Head Diameter : 6.0mm
Outer Diameter : 4.0mm | Inner Diameter : 1.9mm | Pitch : 1.75mm


PRODUCT CODE

LENGTH

MT-SS1740010	10
MT-SS1740012	12
MT-SS1740014	14
MT-SS1740016	16
MT-SS1740018	18
MT-SS1740020	20
MT-SS1740022	22
MT-SS1740024	24
MT-SS1740026	26
MT-SS1740028	28
MT-SS1740030	30
MT-SS1740035	35
MT-SS1740040	40
MT-SS1740045	45
MT-SS1740050	50
MT-SS1740055	55
MT-SS1740060	60
MT-SS1740065	65
MT-SS1740070	70
MT-SS1740075	75
MT-SS1740080	80



KET™ 3.5mm Locking Plate system Instrument Set



KET™ 3.5mm Locking Plate system Instrument Set

Part No:MT-INSLPS002	Description
	Drill Bit Ø2.5mm, Length 125mm, 2 flute, for Quick coupling
	
Part No:MT-INSLPS004	Description
	Drill Bit Ø3.5mm, Length 125mm, 2 flute, for Quick coupling
	
Part No:MT-INSLPS003	Description
	Drill Bit Ø2.8mm, Length 147mm, 2-flute, for Quick coupling
	
Part No:MT-INSLPS044	Description
	Drill Bit Ø2.8mm, Length 200mm, 2-flute, for Quick coupling
	
Part No:MT-INSLPS011	Description
	Double drill Guide 3.5mm & 2.5 mm
	

KET™ 3.5mm Locking Plate system Instrument Set

Part No:MT-INSLPS015

Description

Universal Drill Guide 2.5mm & 3.5mm



Part No:MT-INSLPS023

Description

2.5 Screw Driver shaft with quick coupling for Torque Limiter



Part No:MT-INSLPS022

Description

Manual 1.5Nm Torque Limiter



Part No:MT-INSLPS032

Description

Quick Change Handle for SLPS



Part No:MT-INSLPS024

Description

Screwdriver, Hexagonal, small, A/F 2.5mm, with Groove & Holding Sleeve



KET™ 3.5mm Locking Plate system Instrument Set

Part No:MT-INSLPS028

Description

Depth gauge for Screw 2.7/3.5mm measuring range up to 90mm



Part No:MT-INSLPS031

Description

T-Handle with Quick coupling SLPS



Part No:MT-INSLPS006

Description

Tap for Cortex screw Ø3.5mm Length 115mm for Quick Coupling



Part No:MT-INSLPS012

Description

LPS Drill Sleeve for Drill Bits Ø2.8mm, Length 65mm



Part No:MT-INSLPS020

Description

LPS Drill Sleeve for Drill Bits Ø2.8mm, Length 90mm



KET™ 3.5mm Locking Plate system Instrument Set

Part No:MT-INSLPS026

Description

Wrench-SLPS



Part No:MT-INSLPS017

Description


Guide Wire - Krischner - SLPS



Part No:MT-INSLPS043

Description

Instrument Case with Lid



KET™ Elbow-Clavicle Plate System

2.7mm/3.5mm Superior Anterior Clavicle Plate

Material : SS 316 LVM | Thickness : 3 mm

PRODUCT CODE	HOLES	LENGTH
MT-PS1006094L	6	92
MT-PS1007110L	7	106
MT-PS1008120L	8	120
MT-PS1006094R	6	92
MT-PS1007110R	7	106
MT-PS1008120R	8	120



2.7mm/3.5mm Superior Anterior Clavicle Plate, with Lateral Extension

Material : SS 316 LVM | Thickness : 3 mm

PRODUCT CODE	HOLES	LENGTH
MT-PS1103069L	3	66
MT-PS1104081L	4	80
MT-PS1105094L	5	94
MT-PS1106108L	6	108
MT-PS1107122L	7	122
MT-PS1108136L	8	136
MT-PS1103069R	3	66
MT-PS1104081R	4	80
MT-PS1105094R	5	94
MT-PS1106108R	6	108
MT-PS1107122R	7	122
MT-PS1108136R	8	136



2.7mm/3.5mm Posterolateral Distal Humerus Plate

Material : SS 316 LVM | Thickness : 2.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PS0203055L	3	65
MT-PS0205090L	5	91
MT-PS0207116L	7	117
MT-PS0209143L	9	143
MT-PS0214208L	14	208
MT-PS0203055R	3	65
MT-PS0205090R	5	91
MT-PS0207116R	7	117
MT-PS0209143R	9	143
MT-PS0214208R	14	208



KET™ Elbow-Clavicle Plate System

2.7mm/3.5mm Posterolateral Distal Humerus Plate, with Lateral Support

Material : SS 316 LVM | Thickness : 2.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PS0303055L	3	65
MT-PS0305090L	5	91
MT-PS0307116L	7	117
MT-PS0309143L	9	143
MT-PS0314208L	14	208
MT-PS0303055R	3	65
MT-PS0305090R	5	91
MT-PS0307116R	7	117
MT-PS0309143R	9	143
MT-PS0314208R	14	208



2.7mm/3.5mm Medial Distal Humerus Plate

Material : SS 316 LVM | Thickness : 2.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PS0103058L	3	58
MT-PS0105083L	5	84
MT-PS0107110L	7	110
MT-PS0109149L	9	136
MT-PS0114201L	14	201
MT-PS0103058R	3	58
MT-PS0105083R	5	84
MT-PS0107110R	7	110
MT-PS0109149R	9	136
MT-PS0114201R	14	201



2.7mm/3.5mm LPS Medial Distal Humerus Metaphyseal Plate

Material : SS 316 LVM | Thickness : 2.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PS1507080	7	80
MT-PS1509106	9	106
MT-PS1511132	11	132
MT-PS1513158	13	158
MT-PS1515184	15	184



KET™ Elbow-Clavicle Plate System

LPS Distal Humerus Extra Articular Plate 3.5

Material : SS 316 LVM | Thickness : 2.5 mm

PRODUCT CODE	HOLES	LENGTH
MT-PS4106158L	6	158
MT-PS4108194L	8	194
MT-PS4110230L	10	230
MT-PS4106158R	6	158
MT-PS4108194R	8	194
MT-PS4110230R	10	230



LPS Olecranon Plate 3.5

Material : SS 316 LVM | Thickness : 3.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS0902086L	2	86
MT-PS0904112L	4	112
MT-PS0906138L	6	138
MT-PS0908164L	8	216
MT-PS0902086R	2	86
MT-PS0904112R	4	112
MT-PS0906138R	6	138
MT-PS0908164R	8	164



KET™ 3.5mm Screw

3.5mm Cortex Screw, Self Tapping | Material : SS 316 LVM | Head Diameter : 6.0mm
Outer Diameter : 3.5mm | Inner Diameter : 2.4mm | Pitch : 1.75mm


PRODUCT CODE	LENGTH
MT-SS1035010	10
MT-SS1035012	12
MT-SS1035014	14
MT-SS1035016	16
MT-SS1035018	18
MT-SS1035020	20
MT-SS1035022	22
MT-SS1035024	24
MT-SS1035026	26
MT-SS1035028	28
MT-SS1035030	30
MT-SS1035032	32
MT-SS1035034	34
MT-SS1035036	36
MT-SS1035038	38
MT-SS1035040	40
MT-SS1035042	42
MT-SS1035044	44
MT-SS1035045	45
MT-SS1035046	46
MT-SS1035048	48
MT-SS1035050	50
MT-SS1035055	55
MT-SS1035060	60
MT-SS1035065	65
MT-SS1035070	70
MT-SS1035075	75
MT-SS1035080	80
MT-SS1035085	85
MT-SS1035090	90
MT-SS1035095	95
MT-SS1035100	100
MT-SS1035105	105
MT-SS1035110	110



KET™ 3.5mm Screw

3.5mm Locking Screw, Self Tapping | Material : SS 316 LVM | Head Diameter : 5.0mm
Outer Diameter : 3.5mm | Inner Diameter : 2.9mm | Pitch : 0.8mm

PRODUCT CODE	LENGTH
MT-SS1235010	10
MT-SS1235012	12
MT-SS1235014	14
MT-SS1235016	16
MT-SS1235018	18
MT-SS1235020	20
MT-SS1235022	22
MT-SS1235024	24
MT-SS1235026	26
MT-SS1235028	28
MT-SS1235030	30
MT-SS1235032	32
MT-SS1235034	34
MT-SS1235036	35
MT-SS1235038	36
MT-SS1235040	38
MT-SS1235042	40
MT-SS1235044	42
MT-SS1235046	46
MT-SS1235048	48
MT-SS1235050	50
MT-SS1235052	52
MT-SS1235054	54
MT-SS1235056	56
MT-SS1235058	58
MT-SS1235060	60
MT-SS1235065	65
MT-SS1235070	70
MT-SS1235075	75
MT-SS1235080	80
MT-SS1235085	85



KET™ 3.5mm Screw

4.0 mm Cancellous Screw, Full Thread | Material : SS 316 LVM | Head Diameter : 6.0mm
Outer Diameter : 4.0mm | Inner Diameter : 1.9mm | Pitch : 1.75mm


PRODUCT CODE	LENGTH
MT-SS1640010	10
MT-SS1640012	12
MT-SS1640014	14
MT-SS1640016	16
MT-SS1640018	18
MT-SS1640020	20
MT-SS1640022	22
MT-SS1640024	24
MT-SS1640026	26
MT-SS1640028	28
MT-SS1640030	30
MT-SS1640035	35
MT-SS1640040	40
MT-SS1640045	45
MT-SS1640050	50
MT-SS1640055	55
MT-SS1640060	60
MT-SS1640065	65
MT-SS1640070	70
MT-SS1640075	75
MT-SS1640080	80



KET™ 3.5mm Screw

4.0mm Cancellous Screw, Short Thread | Material : SS 316 LVM | Head Diameter : 6.0mm
Outer Diameter : 4.0mm | Inner Diameter : 1.9mm | Pitch : 1.75mm

PRODUCT CODE	LENGTH
MT-SS1740010	10
MT-SS1740012	12
MT-SS1740014	14
MT-SS1740016	16
MT-SS1740018	18
MT-SS1740020	20
MT-SS1740022	22
MT-SS1740024	24
MT-SS1740026	26
MT-SS1740028	28
MT-SS1740030	30
MT-SS1740035	35
MT-SS1740040	40
MT-SS1740045	45
MT-SS1740050	50
MT-SS1740055	55
MT-SS1740060	60
MT-SS1740065	65
MT-SS1740070	70
MT-SS1740075	75
MT-SS1740080	80



KET™ 3.5mm Screw

Ø 2.7mm Locking Screw, T8 Self Tapping | Material : SS 316 LVM
Head Diameter : 3.5 | Outer Diameter : 2.7mm | Core Diameter : 2.1mm | Pitch : 0.6mm

PRODUCT CODE

LENGTH

MT-SS0727010	10
MT-SS0727012	12
MT-SS0727014	14
MT-SS0727016	16
MT-SS0727018	18
MT-SS0727020	20
MT-SS0727022	22
MT-SS0727024	24
MT-SS0727026	26
MT-SS0727028	28
MT-SS0727030	30
MT-SS0727032	32
MT-SS0727034	34
MT-SS0727036	36
MT-SS0727038	38
MT-SS0727040	40
MT-SS0727042	42
MT-SS0727044	44
MT-SS0727046	46
MT-SS0727048	48
MT-SS0727050	50
MT-SS0727055	55
MT-SS0727060	60



KET™ Elbow and Clavicle system Instrument Set



KET™ Elbow and Clavicle system Instrument Set

Part No:MT-INDR24002

Description

Drill Bit Ø2.0mm, Length 110mm, 2-flute, for Quick Coupling



Part No:MT-INDR24015

Description

Screwdriver shaft, T8, self-holding



Part No:MT-INDR24013

Description

Manual 0.8Nm Torque Limiter



Part No:MT-INDR24020

Description

Quick Change Handle Large with Quick Coupling, Length 138



Part No:MT-INDR24021

Description

Grasping sets of 2.4/2.7 locking screw, Screw Driver T8



Part No:MT-INDR24006

Description

LPS Drill Sleeve 2.7 with Scale up to 30mm, for Drill bits Ø2.0mm



KET™ Elbow and Clavicle system Instrument Set

Part No:MT-INDR24016

Description

Screw Driver Stardrive T8 Self Holding



Part No:MT-INDR24019

Description

Depth Gauge for Screws Ø2.7 to 4.0 mm, measuring range up to 60 mm



Part No:MT-INELCV001

Description

Instrument Case with Lid




KET™ 4.5mm Locking Plate System

4.5mm LPS Narrow Compression Plate

Material : SS 316 LVM | Thickness : 4.2mm


PRODUCT CODE	HOLES	LENGTH
MT-PS2906116	6	116
MT-PS2907134	7	134
MT-PS2908152	8	152
MT-PS2909170	9	170
MT-PS2910188	10	188
MT-PS2911206	11	206
MT-PS2912224	12	224
MT-PS2914260	14	260



4.5mm LPS Broad Compression Plate

Material : SS 316 LVM | Thickness : 5.2mm

PRODUCT CODE	HOLES	LENGTH
MT-PS2606116	6	116
MT-PS2607134	7	134
MT-PS2608152	8	152
MT-PS2609170	9	170
MT-PS2610188	10	188
MT-PS2611206	11	206
MT-PS2612224	12	224
MT-PS2614260	14	260



4.5mm LPS Distal Femur Plate

Material : SS 316 LVM | Thickness : 5.5mm

PRODUCT CODE	HOLES	LENGTH
MT-PS4605156L	5	156
MT-PS4607196L	7	196
MT-PS4609236L	9	236
MT-PS4611276L	11	276
MT-PS4613316L	13	316
MT-PS4605156R	5	156
MT-PS4607196R	7	196
MT-PS4609236R	9	236
MT-PS4611276R	11	276
MT-PS4613316R	13	316



KET™ 4.5mm Locking Plate System

4.5mm LPS Proximal Lateral Tibia Plate

Material : SS 316 LVM | Thickness : 4.5mm


PRODUCT CODE	HOLES	LENGTH
MT-PS4805140L	5	140
MT-PS4807180L	7	180
MT-PS4809220L	9	220
MT-PS4811260L	11	260
MT-PS4813300L	13	300
MT-PS4805140R	5	140
MT-PS4807180R	7	180
MT-PS4809220R	9	220
MT-PS4811260R	11	260
MT-PS4813300R	13	300



4.5mm LPS T Plate

Material : SS 316 LVM | Thickness : 2.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS3004084	4	84
MT-PS3005100	5	100
MT-PS3006116	6	110
MT-PS3007132	7	132
MT-PS3008148	8	148



4.5mm LPS T Buttress Plate

Material : SS 316 LVM | Thickness : 2.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS2804083	4	83
MT-PS2805099	5	99
MT-PS2806115	6	115



KET™ 4.5mm Locking Plate System

4.5mm LPS L Buttress Plate
Material : SS 316 LVM | Thickness : 2.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS2704088L	4	88
MT-PS2705104L	5	104
MT-PS2706120L	6	120
MT-PS2704088R	4	88
MT-PS2705104R	5	104
MT-PS2706120R	6	120



4.5mm LPS Large Metaphyseal Plate
Material : SS 316 LVM | Thickness : 4.5mm

PRODUCT CODE	HOLES	LENGTH
MT-PS0403118	3	118
MT-PS0404136	4	136
MT-PS0405154	5	154
MT-PS0406172	6	172
MT-PS0407190	7	190
MT-PS0408208	8	208
MT-PS0409226	9	226



KET™ 4.5mm Large Fragment Plate System

4.5mm Low Contourable Dual Compression Plate Narrow
Material : SS 316 LVM | Thickness : 4.2mm

PRODUCT CODE	HOLES	LENGTH
MT-PS3702034	2	34
MT-PS3703052	3	52
MT-PS3704070	4	70
MT-PS3705088	5	88
MT-PS3706106	6	106
MT-PS3707124	7	124
MT-PS3708142	8	142
MT-PS3709160	9	160
MT-PS3710178	10	178
MT-PS3711196	11	196
MT-PS3712214	12	214
MT-PS3713232	13	232
MT-PS3714250	14	250
MT-PS3715268	15	268
MT-PS3716286	16	286



4.5mm Low Contourable Dual Compression Plate Broad
Material : SS 316 LVM | Thickness : 5.2mm

PRODUCT CODE	HOLES	LENGTH
MT-PS3606106	6	106
MT-PS3607124	7	124
MT-PS3608142	8	142
MT-PS3609160	9	160
MT-PS3610178	10	178
MT-PS3611196	11	196
MT-PS3612214	12	214
MT-PS3614250	14	250
MT-PS3616286	16	286
MT-PS3618322	18	322



KET™ 4.5mm Screw

5.0mm Locking Screw, Self Tapping | Material : SS 316 LVM | Head Diameter : 7.0mm
Outer Diameter : 5.0mm | Core Diameter : 4.4mm | Pitch : 1.0mm

PRODUCT CODE

LENGTH

MT-SS2850014	14
MT-SS2850016	16
MT-SS2850018	18
MT-SS2850020	20
MT-SS2850022	22
MT-SS2850024	24
MT-SS2850026	26
MT-SS2850028	28
MT-SS2850030	30
MT-SS2850032	32
MT-SS2850034	34
MT-SS2850036	36
MT-SS2850038	38
MT-SS2850040	40
MT-SS2850042	42
MT-SS2850044	44
MT-SS2850046	46
MT-SS2850048	48
MT-SS2850050	50
MT-SS2850055	55
MT-SS2850060	60
MT-SS2850065	65
MT-SS2850070	70
MT-SS2850075	75
MT-SS2850080	80
MT-SS2850085	85
MT-SS2850090	90



5.0mm Locking Cancellous Screw, Self Tapping | Material : SS 316 LVM
Head Diameter : 7.0mm | Outer Diameter : 5.0mm | Core Diameter : 3.2mm | Pitch : 1.0mm

PRODUCT CODE

LENGTH

MT-SS8550055	55
MT-SS8550060	60
MT-SS8550065	65
MT-SS8550075	75
MT-SS8550080	80
MT-SS8550085	85
MT-SS8550090	90




KET™ 4.5mm Screw

4.5mm Cortex Screw, Self Tapping | Material : SS 316 LVM | Head Diameter : 8.0mm
Outer Diameter : 4.5mm | Core Diameter : 3.0mm | Pitch : 1.75mm

PRODUCT CODE

LENGTH

MT-SS2045014	14
MT-SS2045016	16
MT-SS2045018	18
MT-SS2045020	20
MT-SS2045022	22
MT-SS2045024	24
MT-SS2045026	26
MT-SS2045028	28
MT-SS2045030	30
MT-SS2045032	32
MT-SS2045034	34
MT-SS2045036	36
MT-SS2045038	38
MT-SS2045040	40
MT-SS2045042	42
MT-SS2045044	44
MT-SS2045046	46
MT-SS2045048	48
MT-SS2045050	50
MT-SS2045052	52
MT-SS2045054	54
MT-SS2045056	56
MT-SS2045058	58
MT-SS2045060	60
MT-SS2045062	62
MT-SS2045064	64
MT-SS2045066	66
MT-SS2045068	68
MT-SS2045070	70
MT-SS2045072	72
MT-SS2045076	76
MT-SS2045080	80



KET™ 4.5mm Screw

6.5mm Cancellous Screw, 16mm Thread Length | Material : SS 316 LVM
Head Diameter : 8.0mm | Outer Diameter : 6.5mm | Core Diameter : 3.0mm | Pitch : 2.75mm

PRODUCT CODE

LENGTH

MT-SS3465030	30
MT-SS3465035	35
MT-SS3465040	40
MT-SS3465045	45
MT-SS3465050	50
MT-SS3465055	55
MT-SS3465060	60
MT-SS3465065	65
MT-SS3465070	70
MT-SS3465075	75
MT-SS3465080	80
MT-SS3465085	85
MT-SS3465090	90
MT-SS3465095	95
MT-SS3465100	100
MT-SS3465105	105
MT-SS3465110	110



6.5mm Cancellous Screw, 32mm Thread Length | Material : SS 316 LVM
Head Diameter : 8.0mm | Outer Diameter : 6.5mm | Core Diameter : 3.0mm | Pitch : 2.75mm

PRODUCT CODE

LENGTH

MT-SS3565045	45
MT-SS3565050	50
MT-SS3565055	55
MT-SS3565060	60
MT-SS3565065	65
MT-SS3565070	70
MT-SS3565075	75
MT-SS3565080	80
MT-SS3565085	85
MT-SS3565090	90
MT-SS3565095	95




KET™ 4.5mm Screw

6.5mm Cancellous Screw, Full Thread Length | Material : SS 316 LVM
Head Diameter : 8.0mm | Outer Diameter : 6.5mm | Core Diameter : 3.0mm | Pitch : 2.75mm


PRODUCT CODE

LENGTH

MT-SS3365020	20
MT-SS3365025	25
MT-SS3365030	30
MT-SS3365035	35
MT-SS3365040	40
MT-SS3365045	45
MT-SS3365050	50
MT-SS3365055	55
MT-SS3365060	60
MT-SS3365065	65
MT-SS3365070	70
MT-SS3365075	75
MT-SS3365080	80
MT-SS3365085	85
MT-SS3365090	90
MT-SS3365095	95
MT-SS3365100	100



KET™ 4.5mm Locking Plate system Instrument Set



KET™ 4.5mm Locking Plate system Instrument Set

Part No:MT-INLLPS001

Description

Drill bit Ø3.2mm, Length 145mm, 2-flute, for Quick Coupling



Part No:MT-INLLPS003

Description

Drill bit Ø4.5mm, Length 145mm, 2-flute, for Quick Coupling



Part No:MT-INLLPS002

Description

Drill bit Ø4.3mm, Length 270 mm, 2-flute, for Quick Coupling



Part No:MT-INLLPS012

Description

Double Drill Guide 4.5mm & 3.2mm



Part No:MT-INLLPS011

Description


Universal Drill Guide 4.5mm & 3.2mm



KET™ 4.5mm Locking Plate system Instrument Set

Part No:MT-INLLPS025	Description
	Depth Gauge for Screws 4.5/6.5mm,measuring range up to 100mm
	
Part No:MT-INLLPS019	Description
	Torque Limiter Screw Driver 3.5mm Hex Self Holding 4Nm
	
Part No:MT-INLLPS024	Description
	Screwdriver, Hexagonal, A/F 3.5mm, with Groove & Holding Sleeve
	
Part No:MT-INLLPS010	Description
	LPS Drill Sleeve 5.0, for Drill Bits Ø4.3mm
	
Part No:MT-INLLPS017	Description
	Drill Sleeve OD4.3mm & ID2.0mm
	
Part No:MT-INLLPS015	Description
	Ø3.2 Insert Drill Sleeve
	

KET™ 4.5mm Locking Plate system Instrument Set

Part No:MT-INSLPS031	Description
	T-Handle with Quick coupling SLPS
	
Part No:MT-INLLPS005	Description
	Tap for Cortex Screws Ø4.5mm, Length 130 mm for quick coupling
	
Part No:MT-INSLPS026	Description
	Wrench-SLPS
	
Part No:MT-INSLPS017	Description
	Guide Wire - Krischner - SLPS
	
Part No:MT-INLLPS041	Description
	Instrument Case with Lid
	

KET™ DHF & DSP Plates System

130° Dynamic Hip Fixation System, Long Barrel Plate

Material : SS 316 LVM | Thickness : 6.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS5402046	2	46
MT-PS5404078	4	78
MT-PS5405094	5	94
MT-PS5406110	6	110
MT-PS5408142	8	142
MT-PS5410174	10	174
MT-PS5412206	12	206
MT-PS5414238	14	238
MT-PS5416270	16	270



135° Dynamic Hip Fixation System, Long Barrel Plate

Material : SS 316 LVM | Thickness : 6.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS5502046	2	46
MT-PS5504078	4	78
MT-PS5505094	5	94
MT-PS5506110	6	110
MT-PS5508142	8	142
MT-PS5510174	10	174
MT-PS5512206	12	206
MT-PS5514238	14	238
MT-PS5516270	16	270



KET™ DHF & DSP Plates System

140° Dynamic Hip Fixation System, Long Barrel Plate

Material : SS 316 LVM | Thickness : 6.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS5802046	2	46
MT-PS5804078	4	78
MT-PS5805094	5	94
MT-PS5806110	6	110
MT-PS5808142	8	142
MT-PS5810174	10	174
MT-PS5812206	12	206
MT-PS5814238	14	238
MT-PS5816270	16	270



145° Dynamic Hip Fixation System, Long Barrel Plate

Material : SS 316 LVM | Thickness : 6.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS5902046	2	46
MT-PS5904078	4	78
MT-PS5905094	5	94
MT-PS5906110	6	110
MT-PS5908142	8	142
MT-PS5910174	10	174
MT-PS5912206	12	206
MT-PS5914238	14	238



KET™ DHF & DSP Plates System

150° Dynamic Hip Fixation System, Long Barrel Plate

Material : SS 316 LVM | Thickness : 6.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS5704078	4	78
MT-PS5705094	5	94
MT-PS5706110	6	110
MT-PS5708142	8	142
MT-PS5710174	10	174



130° Dynamic Hip Fixation System, Short Barrel Plate

Material : SS | Thickness : 6.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS6202046	2	46
MT-PS6204078	4	78
MT-PS6205094	5	94
MT-PS6206110	6	110
MT-PS6208142	8	142
MT-PS6210174	10	174
MT-PS6212206	12	206
MT-PS6214238	14	238
MT-PS6216270	16	270



KET™ DHF & DSP Plates System

135° Dynamic Hip Fixation System, Short Barrel Plate

Material : SS 316 LVM | Thickness : 6.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS5602046	2	46
MT-PS5604078	4	78
MT-PS5605094	5	94
MT-PS5606110	6	110
MT-PS5608142	8	142
MT-PS5610174	10	174
MT-PS5612206	12	206
MT-PS5614238	14	238
MT-PS5616270	16	270



140° Dynamic Hip Fixation System, Short Barrel Plate

Material : SS 316 LVM | Thickness : 6.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS6004078	4	78
MT-PS6005094	5	94
MT-PS6006110	6	110
MT-PS6008142	8	142
MT-PS6010174	10	174
MT-PS6012206	12	206
MT-PS6014238	14	238
MT-PS6016270	16	270



KET™ DHF & DSP Plates System

145° Dynamic Hip Fixation System, Short Barrel Plate

Material : SS 316 LVM | Thickness : 6.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS6104078	4	78
MT-PS6105094	5	94
MT-PS6106110	6	110
MT-PS6108142	8	142
MT-PS6110174	10	174
MT-PS6112206	12	206
MT-PS6114238	14	238
MT-PS6116270	16	270



150° Dynamic Hip Fixation System, Short Barrel Plate

Material : SS 316 LVM | Thickness : 6.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS6302046	2	46
MT-PS6304078	4	78
MT-PS6305094	5	94
MT-PS6306110	6	110
MT-PS6308142	8	142
MT-PS6310174	10	174
MT-PS6312206	12	206
MT-PS6314238	14	238
MT-PS6316270	16	270



KET™ DHF & DSP Plates System

95° Dynamic Supracondylar Plate

Material : SS 316 LVM | Thickness : 6.0mm

PRODUCT CODE	HOLES	LENGTH
MT-PS3906114	6	114
MT-PS3908146	8	146
MT-PS3910178	10	178
MT-PS3912210	12	210
MT-PS3914242	14	242
MT-PS3916274	16	274



Trochanteric Stabilizing Plate

Material : SS 316 LVM | Thickness : 1.6mm


PRODUCT CODE	HOLES	LENGTH
MT-PS4005138	5	138
MT-PS4006148	6	148



KET™ DHF & DSP Screw

Ø 6.5mm Cancellouse Screw, Full Thread | Material : SS 316 LVM
Head Diameter : 8.0mm | Outer Diameter : 6.5mm | Core Diameter : 3.0mm | Pitch : 2.75mm

PRODUCT CODE	LENGTH
MT-SS3365020	20
MT-SS3365025	25
MT-SS3365030	30
MT-SS3365035	35
MT-SS3365040	40
MT-SS3365045	45
MT-SS3365050	50
MT-SS3365055	55
MT-SS3365060	60
MT-SS3365065	65
MT-SS3365070	70
MT-SS3365075	75
MT-SS3365080	80
MT-SS3365085	85
MT-SS3365090	90
MT-SS3365095	95



KET™ DHF & DSP Screw

Ø 6.5mm Cancellous Screw, 16mm Thread Length | Material : SS 316 LVM
Head Diameter : 8.0mm | Outer Diameter : 6.5mm | Core Diameter : 3.0mm | Pitch : 2.75mm

PRODUCT CODE	LENGTH
MT-SS3465030	30
MT-SS3465035	35
MT-SS3465040	40
MT-SS3465045	45
MT-SS3465050	50
MT-SS3465055	55
MT-SS3465060	60
MT-SS3465065	65
MT-SS3465070	70
MT-SS3465075	75
MT-SS3465080	80
MT-SS3465085	85
MT-SS3465090	90
MT-SS3465095	95



KET™ DHF & DSP Screw

Ø 6.5mm Cancellous Screw, 32mm Thread Length | Material : SS 316 LVM
Head Diameter : 8.0mm | Outer Diameter : 6.5mm | Core Diameter : 3.0mm | Pitch : 2.75mm

PRODUCT CODE

LENGTH

MT-SS3565045	45
MT-SS3565050	50
MT-SS3565055	55
MT-SS3565060	60
MT-SS3565065	65
MT-SS3565070	70
MT-SS3565075	75
MT-SS3565080	80
MT-SS3565085	85
MT-SS3565090	90
MT-SS3565095	95



KET™ DHF & DSP Screw

12.5mm DHF/DSP Lag Screws | Material : SS 316 LVM
Outer Diameter : 12.0mm | Core Diameter : 8.2mm | Pitch : 3.0mm

PRODUCT CODE

LENGTH

MT-SS4022050	50
MT-SS4022055	55
MT-SS4022060	60
MT-SS4022065	65
MT-SS4022070	70
MT-SS4022075	75
MT-SS4022080	80
MT-SS4022085	85
MT-SS4022090	90
MT-SS4022095	95
MT-SS4022100	100
MT-SS4022105	105
MT-SS4022110	110
MT-SS4022115	115
MT-SS4022120	120
MT-SS4022125	125
MT-SS4022130	130
MT-SS4022135	135
MT-SS4022140	140
MT-SS4022145	145




Compression Screw for DHF/DSP | Material : SS 316 LVM
Head Diameter : 8.0mm | Outer Diameter : 4.0mm | Core Diameter : 3.3mm | Pitch : 0.7mm

PRODUCT CODE

LENGTH

MT-SS4140036	36
--------------	----




KET™ DHF & DSP Screw

4.5mm Cortex Screw, Self Tapping | Material : SS 316 LVM | Diameter : 8.0mm |
Outer Diameter : 4.5mm | Inner Diameter : 3.0mm | Pitch : 1.75mm


PRODUCT CODE

LENGTH

MT-SS2045014	14
MT-SS2045016	16
MT-SS2045018	18
MT-SS2045020	20
MT-SS2045022	22
MT-SS2045024	24
MT-SS2045026	26
MT-SS2045028	28
MT-SS2045030	30
MT-SS2045032	32
MT-SS2045034	34
MT-SS2045036	36
MT-SS2045038	38
MT-SS2045040	40
MT-SS2045042	42
MT-SS2045044	44
MT-SS2045046	46
MT-SS2045048	48
MT-SS2045050	50
MT-SS2045052	52
MT-SS2045054	54
MT-SS2045056	56
MT-SS2045058	58
MT-SS2045060	60
MT-SS2045062	62
MT-SS2045064	64
MT-SS2045066	66
MT-SS2045068	68
MT-SS2045070	70
MT-SS2045072	72
MT-SS2045074	74
MT-SS2045076	76
MT-SS2045080	80
MT-SS2045085	85
MT-SS2045090	90
MT-SS2045095	95
MT-SS2045100	100
MT-SS2045105	105
MT-SS2045110	110
MT-SS2045115	115
MT-SS2045120	120
MT-SS2045125	125
MT-SS2045130	130
MT-SS2045135	135
MT-SS2045140	140




KET™ DHF & DSP Plate System Instrument Set



KET™ Dynamic Fixation Plate system Instrument Set

Part No:MT-INLLPS001	Description
	Drill Bit Ø3.2mm, Length 145mm, 2 flute, for Quick coupling
	
Part No:MT-INLLPS003	Description
	Drill Bit Ø4.5mm, Length 145mm, 2 flute, for Quick coupling
	
Part No:MT-INDHFS003	Description
	Drill Bit Ø8.0mm, Length 245mm, 2 flute, for Quick coupling
	
Part No:MT-INDHFS026	Description
	DHF Triple Reamer, Long Reamer
	
Part No:MT-INDHFS027	Description
	DSP Triple Reamer
	
Part No:MT-INDHFS004	Description
	DHF/DSP Tap
	

KET™ Dynamic Fixation Plate system Instrument Set

Part No:MT-INDHFS006	Description
	Tap HB 6.5
	
Part No:MT-INDHFS024	Description
	Simple Drill chuck
	
Part No:MT-INDHFS011	Description
	Double Drill Guide 3.2mm & 6.5mm
	
Part No:MT-INDHFS012	Description
	Double Drill Guide 3.2mm & 3.2mm, for Neutral and Load position
	
Part No:MT-INLLPS012	Description
	Double Drill Guide 3.2mm & 4.5mm
	

KET™ Dynamic Fixation Plate system Instrument Set

Part No:MT-INDHFS014	Description
	DSP Angled Guide 95°
	
Part No:MT-INDHFS008	Description
	DHF Angled Guide 135°
	
Part No:MT-INDHFS009	Description
	DHF/DSP Centering Sleeve
	
Part No:MT-INDHFS007	Description
	DHF/DSP Guide Wire Ø2.5mm with threaded tip with trocar,Length 230mm
	

KET™ Dynamic Fixation Plate system Instrument Set

Part No:MT-INLLPS024

Description

Screwdriver, Hexagonal, A/F 3.5mm, with Groove & Holding Sleeve



Part No:MT-INSLPS031

Description

T-Handle with Quick coupling SLPS



Part No:MT-INSLPS033

Description

Screw Holding Forceps - SLPS



Part No:MT-INLLPS025

Description

Depth Gauge for Screws 4.5/6.5mm, measuring range up to 100mm



Part No:MT-INDHFS020

Description

DHF/DSP Direct Measuring Device



KET™ Dynamic Fixation Plate system Instrument Set

Part No:MT-INDHFS028

Description

DHF/DSP Impactor



Part No:MT-INDHFS015

Description

DHF Wrench for one-step insertion Length 230mm



Part No:MT-INDHFS019

Description

DHF Connecting Screw, for DHF Wrench



Part No:MT-INDHFS017

Description

DSP Wrench for one-step insertion Length 239mm



KET™ Dynamic Fixation Plate system Instrument Set

Part No:MT-INDHFS018

Description

DSP Connecting Screw, for DSP Wrench



Part No:MT-INDHFS010

Description

Cleaning guide pin Ø2.8



Part No:MT-INDHFS029

Description

DHF/DSP Instrument Case

