


One Step qualitative immunochromatographic assay for detection of antibodies against HIV - 1 and HIV - 2 in human serum, plasma and whole blood (Fingerstick & Venous)




SPECIFICATION

Type	Rapid immunochromatographic test kit
Detected Species	IgM, IgG and IgA Antibodies
Specimen	Serum, Plasma, Whole blood (Fingerstick & Venous)
Specimen Volume	10 µl Serum/ Plasma or 20 µl Whole Blood
Buffer Volume	3 drops
Specimen transfer device	Capillary Tube
Read Time	20 minutes
Storage Temperature	2-30°C
Shelf Life	24 months
Diagnostic Sensitivity	100 %
Diagnostic Specificity	100 %

TEST PROTOCOL



INTERPRETATION



ORDERING INFORMATION

Product Code	Material Description	Pack Size
HVWRPD-01	MERISCREEN HIV 1-2 WB	30 Test
HVWRPD-02	MERISCREEN HIV 1-2 WB	40 Test
HVWRPD-06	MERISCREEN HIV 1-2 WB	50 Test
HVWRPD-07	MERISCREEN HIV 1-2 WB	10 Test
HVWRPD-08	MERISCREEN HIV 1-2 WB	100 Test
HVWRPD-09	MERISCREEN HIV 1-2 WB	25Test
HVWRPD-10	MERISCREEN HIV 1-2 WB (Auto Safety Lancet)	30 Test
HVWRPD-11	MERISCREEN HIV 1-2 WB (Auto Safety Lancet)	60 Test
HVWRPD-12	MERISCREEN HIV 1-2 WB (Auto Safety Lancet)	40 Test



Meril Diagnostics Pvt. Ltd.

D1- D3. Second Floor. Meril Park. Survey No. 135/2/B & 174/2, Muktanand Marg, Chala, Vapi - 396191, India

+91 9624025718 | Email: mahendra.s@merillife.com | askinfo@merillife.com