


MERISCREEN Malaria Pf HRP-II Ag



SPECIFICATION


Type	Rapid immunochromatographic test kit
Detected Species	<i>P. falciparum</i>
Specimen	Whole Blood
Specimen transfer device	Inverted Cup
Read Time	20 minutes
Storage Temperature	1-40°C
Shelf Life	24 months
Limit of Detection (LOD)	≥50 parasites/μl

MERISCREEN Malaria Pf HRP-II Ag


Meril

Diagnostics

TEST PROTOCOL



INTERPRETATION



ORDERING INFORMATION

PRODUCT CODE	DESCRIPTION	Pack Size
RPWMPH-02	MERISCREEN Malaria Pf HRP-II Ag	10 Tests
RPWMPH-03	MERISCREEN Malaria Pf HRP-II Ag	25 Tests
RPWMPH-04	MERISCREEN Malaria Pf HRP-II Ag	50 Tests
RPWMPH-05	MERISCREEN Malaria Pf HRP-II Ag (POCT)	1 Test/kit
RPWMPH-06	MERISCREEN Malaria Pf HRP-II Ag (POCT)	1 Test/kit x 10 each
RPWMPH-07	MERISCREEN Malaria Pf HRP-II Ag (POCT)	1 Test/kit x 25 each

WBRO/IMG/RPWMPH02/01 DATE:16/12/2024



Meril Diagnostics Pvt. Ltd.

D1- D3, Second Floor, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg, Chala, Vapi - 396191, India

+91 9624025718 | Email: mahendra.s@merillife.com | askinfo@merillife.com