

MOZEC™ PEB PTA

Paclitaxel Eluting PTA Balloon Catheter

Meril

More to Life

Extend your reach (of drug treatment)
Longest length 120mm



- Avoid the need for multiple balloons to cover longer lesion
- Lesser exposure time and faster procedure outcome

Securing Outcomes

- Proprietary drug adhesion technology ensures no drug loss during loading due to physical content
- Prospects the integrity of the drug and excipient

Enhanced Navigation: Mozec PTA platform

- Designed on Mozec PTA balloon - a very navigable and well known balloon platform
- PEB available in 0.035" OTW system



MOZEC™ PEB TECHNICAL SPECIFICATIONS

Catheter System	: Over-the-wire (OTW)
Drug Type :	: Paclitaxel
Drug Dosage	: 3 µg/mm ²
Usable Catheter Length	: 80 cm & 135 cm
Max. Guidewire Compatibility	: 0.035"
Balloon Diameter (mm)	: 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00 & 10.00 mm
Balloon Length (mm)	: 30, 40, 50, 60, 80 & 120
Min. Sheath Compatibility	: 5 Fr to 8 Fr (Refer product label/IFU for more details)
Balloon RO Markers	: 2 - Platinum/Iridium
Shaft Diameter	: 5 Fr
Nominal Pressure(NP)	: 6-7 atm depending upon diameter and length of balloon(Refer product label/IFU for more details)
Rated Burst Pressure (RBP)	: 11-14 atm depending upon diameter and length of balloon(Refer product label/IFU for more details)

Mozec™ PEB PTA - ORDERING INFORMATION

Usable Catheter length 80 cm						
Balloon Diameter (mm)	Balloon Length (mm)					
	30	40	50	60	80	120
3.00	MO35P030030A	MO35P030040A	MO35P030050A	MO35P030060A	MO35P030080A	MO35P030120A
4.00	MO35P040030A	MO35P040040A	MO35P040050A	MO35P040060A	MO35P040080A	MO35P040120A
5.00	MO35P050030A	MO35P050040A	MO35P050050A	MO35P050060A	MO35P050080A	MO35P050120A
6.00	MO35P060030A	MO35P060040A	MO35P060050A	MO35P060060A	MO35P060080A	MO35P060120A
7.00	MO35P070030A	MO35P070040A	MO35P070050A	MO35P070060A	MO35P070080A	MO35P070120A
8.00	MO35P080030A	MO35P080040A	MO35P080050A	MO35P080060A	MO35P080080A	MO35P080120A
9.00	MO35P090030A	MO35P090040A	MO35P090050A	MO35P090060A	MO35P090080A	-
10.00	MO35P100030A	MO35P100040A	MO35P100050A	MO35P100060A	MO35P100080A	-

Usable Catheter length 135 cm						
Balloon Diameter (mm)	Balloon Length (mm)					
	30	40	50	60	80	120
3.00	MO35P030030B	MO35P030040B	MO35P030050B	MO35P030060B	MO35P030080B	MO35P030120B
4.00	MO35P040030B	MO35P040040B	MO35P040050B	MO35P040060B	MO35P040080B	MO35P040120B
5.00	MO35P050030B	MO35P050040B	MO35P050050B	MO35P050060B	MO35P050080B	MO35P050120B
6.00	MO35P060030B	MO35P060040B	MO35P060050B	MO35P060060B	MO35P060080B	MO35P060120B
7.00	MO35P070030B	MO35P070040B	MO35P070050B	MO35P070060B	MO35P070080B	MO35P070120B
8.00	MO35P080030B	MO35P080040B	MO35P080050B	MO35P080060B	MO35P080080B	MO35P080120B
9.00	MO35P090030B	MO35P090040B	MO35P090050B	MO35P090060B	MO35P090080B	-
10.00	MO35P100030B	MO35P100040B	MO35P100050B	MO35P100060B	MO35P100080B	-

Sheath Compatibility(OTW 0.035")

5Fr for balloon diameters 3.00, 4.00 and 5.00 mm
6Fr for balloon diameters 6.00 and 7.00 mm
7Fr for balloon diameter 8.00 mm
8Fr for balloon diameters 9.00 and 10.00 mm

CE 1783

Contact your country local Meril sales representative for availability of sizes.
Mozec™ PEB PTA is a registered trademark of Meril Life Sciences Pvt. Ltd.
Mozec™ PEB PTA is not approved and not available for sale in USA.