



Meril

MY♡SIST™

Intra Aortic Balloon Catheter Kit (8F)



For
**Cardiogenic Shock, CABG, Hemodynamically
Unstable Infarction, CVTS*...**

MYOSIST™

Intra Aortic Balloon Catheter Kit (8F)

Key Features:

♥ Highly Durable & Flexible Balloon Membrane

Tested for 1 million inflation deflation cycle for functionality and performance integrity.

♥ Radio-opaque Catheter Tip & Markers

Easy to locate the Myosist™ deployment

♥ Optimized Insertion Length

With 710mm insertion length


♥ Inner Lumen Diameter

0.028 inch diameter provides a medium to measure excellent arterial pressure


Exceptional Benefits:


- ♥ Greater hemodynamic support for every patient
- ♥ More blood volume displacement
- ♥ More diastolic augmentation
- ♥ More systolic unloading

MYOSIST™ Accessories Tray



MYOSIST™ Catheter Tray





Balloon Membrane Sizing

Balloon Volume (cc)	Balloon Dimensions		Patient Height (cm)
	Diameter (mm)	Length (mm)	
30	14.5	225	145-155
35	15	227	155-165
40	15	259	> 165

Ordering Information

MYOSIST IAB Catheter kit 8F



Technical Specifications

Balloon Volume (CC)	30	35	40
Balloon Diameter(mm)	14.5	15	15
Balloon Length(mm)	225	227	259
Balloon Catheter Length(mm)		870	
Balloon Catheter Usable Length(mm)		710	
Catheter Profile(French)		8	
Inner Lumen ID (inch)		0.028	
Guidewire Compatibility (inch)		0.025	
Nominal Pressure (PSI)		2	
Maximum Dilation Pressure (PSI)		5	
Console Compatibility	Datascope : System 97, 98, 98XT, CS 100, CS 300, CardioSave Arrow : ACAT, AutoCAT, AutoCAT 2, AutoCAT 2 Wave, KAAT II		
Balloon Material	Polyurethane		

Myosist™ is a registered trademark of Meril Life Sciences Pvt. Ltd.



More to Life

Manufacturer:

Meril Life Sciences Pvt. Ltd.

Survey No. 135/139, Bilakhia House,
Muktanand Marg, Chala, Vapi - 396 191.
Gujarat. India.
T +91 260 240 8000

Subsidiary companies:

Meril Life Sciences Pvt. Ltd.

301, A-Wing, Business Square,
Chakala, Andheri Kurla Road,
Andheri East, Mumbai 400 093
T +91 22 39350700
F +91 22 39350777

Meril, Inc.

2436 Emrick Boulevard,
Bethlehem, PA 18020
T +610 500 2080
F +610 317 1672

Meril South America

Doc Med LTDA
Al. dos Tupiniquins,
1079 - Cep: 04077-003 - Moema.
Sao Paulo, Brazil.
T +55 11 3624 5935
F +55 11 3624 5936

Meril GmbH.

Bornheimer Strasse 135-137,
D-53119 Bonn.
Germany.
T +49 228 7100 4000
F +49 228 7100 4001

Meril Tibbi Cihazlar

Meril Tibbi Cihazlar İmalat ve Ticaret A.Ş.
İçerenköy Mah. Çetinkaya Sok.
Prestij Plaza No:28
Kat:4 Ataşehir, 34752
İstanbul / Turkey
T +90 216 641 44 24
F +90 216 641 44 25

Meril China Co. Ltd.

2301b,23f, Lixin Plaza,
no 90, South Hubin Road,
Xiamen, China
T 0086-592-5368505
F 0086-592-5368519

Meril SA Pty. Ltd.

102, 104, S101 and S102,
Boulevard West Office Park,
142 Western Service Road,
Erf 813 Woodmead Extension 17
Sandton, Johannesburg – 2191
South Africa
T +27 11 465-2049
F +27 86 471 7941

EU representative.

Obelis S.A.
Bd, General Wahis 53,
1030, Brussels, Belgium.
T +32 2 732 5954
F +32 2 732 6003
E mail@obelis.net

Meril Medical LLC.

Nauchnyi Proezd 19,
Moscow, Russia – 117 246.
Office - +7 495 772 7643

E askinfo@merillife.com
W www.merillife.com