

Meril

MEGASORBTM

Absorbable Polyglycolic Acid Surgical Suture





**Alternative choice for
Soft tissue approximation**

ADVANTAGES

- High tensile strength
- Excellent knotting ability
- Smooth passage through tissue
- Easy handling

PRODUCT OVERVIEW

Structure	: Multifilament, Braided
Colour	: Dyed, Undyed
	
Chemical composition	: Polyglycolic acid
Coating	: Polycaprolactone and Calcium Stearate
Sizes	: USP 8-0 (0.4 metric) to USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilisation	: Ethylene oxide



Meril Endo Surgery Pvt.Ltd.

Meril Park, F3 (E1-E3), Survey No. 1575 and Meril Park 2,
F0, F1, F2 & F4, Survey No. 879, Muktanand Marg,
Chala, Vapi, Valsad, Gujarat, 396191, India