

MeriTe

Cu 250 | Cu 375 | Cu-T 380A
Intrauterine Contraceptive Device

Meril

Endo-Surgery




The Future is Mine

MeriTe | Cu-T 380A

Intrauterine Contraceptive Device

10 years

- Large copper surface area of 380mm²
- The body of the device contains barium sulfate making it visible to X-rays
- Comes with high density polyethylene thread
- Has a simple insertion & removal technique
- More than 99% efficiency




MeriTe | Cu 375

Intrauterine Contraceptive Device

5 years

- U-shaped design for better retention
- Provides freedom from regular check-ups
- Has a simple insertion & removal technique
- Comes with Nylon thread for greater comfort
- Also available in sleek model




MeriTe | Cu 250

Intrauterine Contraceptive Device

3 years

- U-shaped design for better retention
- Provides freedom from regular check-ups
- Has a simple insertion & removal technique
- Comes with Nylon thread for greater comfort
- Also available in sleek model



MeriTe | Cu-T 380A
Renders contraception for 10 years



MeriTe | Cu 375
Renders contraception for 5 years



MeriTe | Cu 250
Renders contraception for 3 years

MERIL ENDO-SURGERY PVT. LTD.

Meril Park, F3 (E1-E3), Survey No. 1575 and
Meril Park 2, F0, F1, F2 & F4, Survey No. 879,
Muktanand Marg, Chala, Vapi, Valsad,
Gujarat, 396191, India

India Sales and Marketing, Mumbai

601-612, Midas, Sahar Plaza Complex, M.V. Road,
J.B Nagar, Andheri (East), Mumbai 400 059, India.
Tel.: +91 22 307 32 400 Email: endosurgery@merillife.com,
Web: www.merillife.com



ME/BRC/IUD/001 Rev.01