

Meril

Endo-Surgery

MERIGLU™

High Viscosity Tissue Adhesive
N-Butyl Cyanoacrylate



MERIGLU™

Topical Skin Adhesive is a waterproof formula consisting of monomeric N-Butyl =Cyanoacrylate, which polymerises quickly when comes in contact with tissue fluid.

MERIGLU™ is intended for topical application to close easily approximated skin edges of wounds from surgical incisions and trauma-induced lacerations.



**WE GOT YOU
GLUED to the
best outcomes.**

MBS™
Meril Biosurgical Solutions

MERIGLUTM

High Viscosity Tissue Adhesive
N-Butyl Cyanoacrylate



MERIGLUTM

Topical Skin Adhesive is a waterproof formula consisting of monomeric N-Butyl Cyanoacrylate, which polymerises quickly when comes in contact with tissue fluid.

MERIGLUTM is intended for topical application to close easily approximated skin edges of wounds from surgical incisions

Advantages

**Faster Set
time <
50 seconds**

**Patient Comfort
and Appreciable
Cosmesis**

Usage Indications

- General surgery – laparoscopy, hernia, appendicitis, vascular surgery
- Plastic surgery – cleft lip, excision of scars, skin graft
- Gynecology – caesarean section, hysterectomy, tubectomy
- Orthopedic – tendon repair, bone graft, cartilage binding
- Surgical oncology – thyroid and breast malignancy surgery
- Pediatrics – management of pediatric lacerations
- Dental surgery
- ENT surgery
- Trauma Care



FAST



RELIABLE



SECURE

**WE GOT YOU
GLUED to the
best outcomes.**