

Meril

Take a leap with

PINIONTM

Knotless Suture

Swift | Facile | Assured
suturing experience...

ME/BRC/024 Rev.02

PINION Knotless Suture

The Pinion knotless suture is a thorough going extensive mechanism of suturing which terminates the need of tying knots in the suture thereby giving a faster and more secure wound closure.

- Absorbable, monofilament, knotless suture
- Shaped Barbs with Dual-axis cut
- Available in uni-directional and bi-directional
- Available in Less Memory Tray to avoid memory impact on sutures


PINION IS IDEAL FOR

- ✓ Simple closures
- ✓ Multilayered closures
- ✓ Closure of high-tension wounds

In a variety of surgical settings



Assured- The suture loop delivers a steadfast assurance, preventing any possibility of wound dehiscence



Facile- The premium needle facilitates smooth navigation through challenging tissues, ensuring effortless passage with ease

SWIFT- PINION Knotlesse's Suture uniquely placed barbs ensure swift glide through tissues, enabling atraumatic wound closure



The design of the barb's tip is carefully crafted to prevent any trauma to the peripheral tissue, ensuring a gentle and non-disruptive interaction.

Pinon's knotless Suture Dual-Axis cuts prevent tissue trauma and minimize dragging for precise surgical outcomes.

Features & Benefits of Pinion

-  **Uniquely shaped barbs**
Enhanced tissue engagement
-  **Pre-knotted loop**
Time efficiency
-  **High tensile strength**
Reduce risk of wound dehiscence



-  **Knotless suture material**
Minimal tissue trauma
-  **Uni & Bi-directional**
Efficient closure & tissue alignment
-  **Airtight Seal with Barb Sutures**
Enhanced Infection Prevention

As compared to a conventional suturing technique, which takes longer time; Pinion - Knotless suture is proven to give a Swift, Facile and an Assured suturing experience. Hence Pinion- Knotless sutures can be an ideal choice for an effectual wound closure.



Swift | Facile | Assured
suturing experience...

FG Code	Suture Type & Length	Size	Dimension	Needle Description
PNDS010307	BARB(PDO VIOLET)15CM	1	37MM	1/2 CIRCLE ROUND BODY UD
PNDS010201	BARB(PDO VIOLET)45CM	1	36MM	1/2 CIRCLE REVERSE CUTTING UD
PNDS010317	BARB(PDO VIOLET)30CM	1	37MM	1/2 CIRCLE TAPER POINT UD
PNDS100306	BARB(PDO VIOLET)15CM	1-0	37MM	1/2 CIRCLE ROUND BODY UD
PNDS100316	BARB(PDO VIOLET)30CM	1-0	37MM	1/2 CIRCLE ROUND BODY UD
PNDS100326	BARB(PDO VIOLET)45CM	1-0	37MM	1/2 CIRCLE ROUND BODY UD
PNDS100336	BARB(PDO VIOLET)60CM	1-0	37MM	1/2 CIRCLE ROUND BODY UD
PNDS100316	BARB(PDO VIOLET)30CM	1-0	37MM	1/2 CIRCLE ROUND BODY UD
PNDS10411BD	BARB(PDO VIOLET)30X30CM	1-0	36MM	1/2 CIRCLE ROUND BODY BD
PNDS200315	BARB (PDO VIOLET) 30CM	2-0	37MM	1/2 CIRCLE ROUND BODY UD
PNDS200403	BARB(PDO VIOLET)45CM	2-0	37MM	1/2 CIRCLE ROUND BODY UD
PNDS200615	BARB(PDO VIOLET)30CM	2-0	26MM	1/2 CIRCLE ROUND BODY UD
PNDS200408	BARB(PDO VIOLET)45CM	2-0	26MM	1/2 CIRCLE ROUND BODY UD
PNCS300625	BARB(PGA-PCL UND)30CM	2-0	36MM	1/2 CIRCLE ROUND BODY UD
PNCS202105	BARB(PGA-PCL UND)15CM	2-0	37MM	1/2 CIRCLE ROUND BODY UD
PNCS202145	BARB(PGA-PCL UND)25CM	2-0	27MM	1/2 CIRCLE ROUND BODY UD
PNDS202145	BARB(PDO VIOLET)25CM	2-0	27MM	1/2 CIRCLE ROUND BODY UD
PNCS20401BD	BARB(PGA-PCL UND)36X36CM	2-0	36MM	1/2 CIRCLE ROUND BODY BD
PNDS300804	BARB(PDO VIOLET)15CM	3-0	17MM	1/2 CIRCLE ROUND BODY UD
PNDS300814	BARB(PDO VIOLET)30CM	3-0	17MM	1/2 CIRCLE ROUND BODY UD
PNCS300014	BARB(PGA-PCL UND)30CM	3-0	19MM	3/8 CIRCLE REVERSE CUTTING UD
PNCS300024	BARB(PGA-PCL UND)45CM	3-0	19MM	3/8 CIRCLE REVERSE CUTTING UD
PNCS300134	BARB(PGA-PCL UND)58CM	3-0	24MM	3/8 CIRCLE REVERSE CUTTING UD
PNCS300604	BARB(PGA-PCL UND)15CM	3-0	26MM	1/2 CIRCLE ROUND BODY UD
PNCS300614	BARB(PGA-PCL UND)30CM	3-0	26MM	1/2 CIRCLE ROUND BODY UD
PNCS300624	BARB(PGA-PCL UND)45CM	3-0	26MM	1/2 CIRCLE ROUND BODY UD
PNDS300604	BARB(PDO VIOLET)15CM	3-0	26MM	1/2 CIRCLE ROUND BODY UD
PNCS300804	BARB(PGA-PCL UND)15CM	3-0	17MM	1/2 CIRCLE ROUND BODY UD
PNDS300024	BARB(PDO VIOLET)45CM	3-0	19MM	3/8 CIRCLE REVERSE CUTTING UD
PNCS30408BD	BARB(PGA-PCL UND)30X30CM	3-0	19MM	1/2 CIRCLE ROUND BODY BD
PNDS400803	BARB(PDO VIOLET)15CM	4-0	17MM	1/2 CIRCLE ROUND BODY UD
PNCS400033	BARB(PGA-PCL UND)58CM	4-0	19MM	3/8 CIRCLE REVERSE CUTTING UD
PNCS40413BD	BARB(PGA-PCL UND)30X30CM	4-0	26MM	1/2 CIRCLE ROUND BODY BD
PNDS010201	BARB(PDO VIO)1 x 45-36MM HC RC UD	1	36MM	1/2 CIRCLE REVERSE CUTTING UD
PNDS010201XL	BARB(PDO VIO)1 x 60-36MM HC RC UD	1	36MM	1/2 CIRCLE REVERSE CUTTING UD



Endo surgery Pvt. Ltd
Meril Park, F3 (E1-E3),Survey No. 1575 and
Meril Park 2, F0, F1, F2 & F4, Survey No. 879,
Muktanand Marg, Chala,Vapi, Valsad,
Gujarat, 396191, India

India Sales & Marketing
Meril Endo-Surgery Pvt. Ltd. 601, Midas,
Sahar Plaza J.B. Nagar, Andheri East, Mumbai 400 059
Maharashtra, India. T : +91 22 307 32400