

# Product Specifications

Intended Use	3 dot Rapid test device with flow through technology for detection of Anti HCV
Principle	Sandwich Immunoassay
Technique	Indirect Immunoassay
Detection	Flow Through
Capture	Antibodies to HCV
Conjugate	Highly purified Recombinant antigens - Core, NS3, NS4 and NS5
Type of Specimen required	Serum or Plasma
Specimen Volume	50 ml (1 Drop)
Result interpretation time	3 minutes
Diagnostic Sensitivity	>99%
Diagnostic Specificity	> 98%
Evaluations	National Institute of Biologicals

## Size Chart

Ordering Information		
PRODUCT CODE	MATERIAL DESCRIPTION	PACK SIZE
HCDRPD-01	Tredro Ab	10 Tests
HCDRPD-03	Tredro Ab	50 Tests