

Product Specifications

Intended Use	3 Dot rapid test device for detection of antibodies against HIV -1 including subtype C and HIV -2 in human serum and plasma
Principle	Indirect Immunoassay
Technique	Flow Through
Detection	Antibodies to HIV -1 and HIV -2
Capture	HIV 1 Region - Recombinant and Synthetic peptides of highly immunodominant region of gp 41 and gp 120 HIV 2 Region - Recombinant and Synthetic peptides of highly immunodominant region of gp 36
Conjugate	Protein A gold Conjugate specific for Fc portion
Type of Specimen required	Serum or Plasma
Specimen Volume	50 ml (1 Drop)
Result interpretation time	3 minutes
Diagnostic Sensitivity	100%
Diagnostic Specificity	> 99.5%
Evaluations	National Institute of Biologicals; National AIDS Research Institute

Size Chart

Ordering Info		
PRODUCT CODE	MATERIAL DESCRIPTION	PACK SIZE
HVDRPD-01	Tredro HIV 1-2 Ab	10 Tests
HVDRPD-03	Tredro HIV 1-2 Ab	50 Tests