

Eliro | Trocar

Trocar sizes - 5mm, 10mm & 12mm



Primary Seal



Bladeless - Tip



Secondary Seal

- **Blunt Bladeless tip:**


- Low Insertion force with simple twisting motion
- Separates rather than cuts along tissue fibers, pushing tissue and vessels away
- Contributes to faster tissue healing

- **Durable Primary Seal:**

- Makes the first entry safe
- Able to withstand pressure along with a long lasting performance
- Accommodates a wide range of instrument sizes

- **Internal Secondary Seal:**

- Facilitates small specimen retrieval, Hernia mesh introduction and Laparoscopic suturing



Code	Product Description	Type Sleeve
ELR05	Bladeless Trocar, 05mm diameter, 100mm Length	Stability
ELR10	Bladeless Trocar, 10mm diameter, 100mm Length	Stability
ELR12	Bladeless Trocar, 12mm diameter, 100mm Length	Stability

Eliro | Trocar-Kit

Trocar-Kit sizes - 5mm, 10mm & 12mm

- † Bladeless Trocars with Tissue Separating Flap
- † One stop solution for all Laparoscopic cases
- † Common Trocar set for Basic Cases – Lap Chole, Appendectomy, TLH, Hernia repair etc
- † Less of Inventory
- † Economical Solution

PACKING



OBTURATOR

SLEEVES/CANNULA

SLEEVES/CANNULA



Code	Description	Type
ELR0505	Two 05mm*100mm SLEEVE/CANNULA and one 05mm OBTURATOR	5MM TROCAR-KIT
ELR1010	Two 10mm*100mm SLEEVE/CANNULA and one 10mm OBTURATOR	10MM TROCAR-KIT
ELR1212	Two 12mm*100mm SLEEVE/CANNULA and one 12mm OBTURATOR	12MM TROCAR-KIT

MERIL ENDO-SURGERY PVT. LTD.

Corporate Headquarters, Vapi

E1-E3, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg,
Chala, Vapi- 396 191, Gujarat, India. Tel.: +91 260 3052 100
Email: endosurgery@merillife.com, Web: www.merillife.com

India Sales & Marketing, Mumbai

601-612, Midas, Sahar Plaza Complex, M.V. Road J.B Nagar,
Andheri (East), Mumbai 400 059, India. Tel.: +91 22 307 32 400
Email: endosurgery@merillife.com,
Web: www.merillife.com