

MIRUSTM

Mirus Endocutter & Reloads



MIRUSTM




Mirus Endocutter & Reloads




Mirus™ Endocutter places two, double or triple-staggered rows of titanium staple (ASTM F 67/ISO 5832-2) and simultaneously divide the tissue between the two rows. The devices are available in 3 different shaft lengths and furnished without a cartridge. Cartridges/Reloads for Endocutter are available in sizes that correspond to the thickness of tissue that is to be transected or resected. The cartridge/reload for Mirus™ Endocutter are supplied separately.

The instruments' safety lock out feature is designed to prevent unintended firing and the retaining cap on the reloads protects the staple leg points during transportation and shipping.

The device with small shaft is used in open surgery whereas the device with medium and large shaft is used in laparoscopic surgery.



Varied Staple Height



New Knife with Every Reload

PRODUCT SPECIFICATIONS

Product Code	Description	Shaft Length
MEC065	MIRUS Endocutter - Short	65 mm
MEC160	MIRUS Endocutter - Medium	160 mm
MEC255	MIRUS Endocutter - Long	255 mm
HMECS45	MIRUS Endocutter 45mm Small	290 mm
HMECM45	MIRUS Endocutter 45mm Medium	340 mm
HMECL45	MIRUS Endocutter 45mm Long	440 mm
HMECS60	MIRUS Endocutter 60mm Small	305 mm
HMECM60	MIRUS Endocutter 60 mm Medium	355 mm
HMECL60	MIRUS Endocutter 60mm Long	455 mm

Code	Trio Reload	Description	Open Staple Height (mm) Inner to Outer
MECR45-2.5	White	MIRUS Endocutter Reload 45mm	2.5
MECR45-3.5	Blue	MIRUS Endocutter Reload 45mm	3.5
MECR45-4.8	Green	MIRUS Endocutter Reload 45mm	4.8
MECR60-2.5	White	MIRUS Endocutter Reload 60mm	2.5
MECR60-3.5	Blue	MIRUS Endocutter Reload 60mm	3.5
MECR60-4.8	Green	MIRUS Endocutter Reload 60mm	4.8
MECRW45	White	MIRUS Endocutter Reload 45mm	2.6, 2.6, 3.0
MECRP45	Purple	MIRUS Endocutter Reload 45mm	3.5, 3.5, 4.0
MECRW60	White	MIRUS Endocutter Reload 60mm	2.6, 2.6, 3.0
MECRP60	Purple	MIRUS Endocutter Reload 60mm	3.5, 3.5, 4.0
MECRB60	Black	MIRUS Endocutter Reload 60mm	4.5, 4.5, 4.8
HMECR45WE	White	MIRUS Endocutter Reload 45mm	2.5
HMECR45BE	Blue	MIRUS Endocutter Reload 45mm	3.5
HMECR45GD	Gold	MIRUS Endocutter Reload 45mm	3.8
HMECR45GE	Green	MIRUS Endocutter Reload 45mm	4.1
HMECR45BK	Black	MIRUS Endocutter Reload 45mm	4.4
HMECR60WE	White	MIRUS Endocutter Reload 60mm	2.5
HMECR60BE	Blue	MIRUS Endocutter Reload 60mm	3.5
HMECR60GD	Gold	MIRUS Endocutter Reload 60mm	3.8
HMECR60GE	Green	MIRUS Endocutter Reload 60mm	4.1
HMECR60BK	Black	MIRUS Endocutter Reload 60mm	4.4

Meril's Global Presence

Meril Endo Surgery Pvt. Ltd.

Third Floor E1-E3, Meril Park, Survey Number
135/2/B & 172/2 Muktanand Marg, Chala,
Vapi 396 191, Gujarat. India.
T : +91 260 240 8000
F : +91 260 305 2125
E : askinfo@merillife.com

Meril Life Sciences Pvt. Ltd

301, A - Wing, Business Square, Chakala,
Andheri Kurla Road, Andheri East,
Mumbai 400 093
T : +91 22 5035 4700

Meril, Inc.

2436 Emrick Boulevard, Bethlehem, PA - 18020
T : (610)500-2080
F : (610)317-1672

Meril South America

Doc Med LTDA, Al. Alameda dos
Maracatins 992, Edificio Aquarius, Andar 8,
Bloco B, Moema - CEP -04089001,
Sao Paulo, SP, Brazil.
T : +55 11 2506 3624
F : +55 11 3624 5936

Meril GmbH.

Bornheimer Strasse 135-137, D-53119
Bonn. Germany.
T : +49 228 7100 4000
F : +49 228 7100 4001

Meril Tibbi Cihazlar

İmalat ve Ticaret A.Ş. İçerenköy Mah.
Çetinkaya Sok. Prestij Plaza No:28,
Kat:4 Ataşehir, 34752 İstanbul / Turkey
T : +90 216 641 44 24
F : +90 216 641 44 25

Meril China Co. Ltd.

2301b,23f, Lixin Plaza,no 90,
South Hubin Road, Xiamen, China.
T : 0086-592-5368505
F : 0086-592-5368519

Meril SA Pty Ltd

102, 104, S101 and S102, Boulevard
West Office Park, 142 Western Service Road,
Erf 813 Woodmead Extension 17 Sandton,
Johannesburg – 2191 South Africa
T : +27 11 465-2049
F : +27 86 471 7941

Meril Medical LLC.

Nauchnyi Proezd 19, Moscow ,
Russia – 117 246.
Office - +7 495 772 7643

EU Representative.

Obelis S.A.Bd, General Wahis 53,1030,
Brussels, Belgium.
T : +32 2 732 5954
F : +32 2 732 6003
E : mail@obelis.net

Meril Bangladesh Pvt. Ltd.

Union Heights-01, 55-2, Bir Uttam Qazi
Nuruzzaman Sarak, West Panthapath,
Dhaka-1205, Bangladesh
T: + 88 02 9122714/24

Meril Australia Pvt. Ltd.

Level 14, 60 Martin Place Sydney
NSW 2000, Australia
T : +61 28248 5800
F : +61 28248 5899

Meril UK

2nd Floor, 55 Ludgate Hill,
London UK, EC4M 7JW