



MIRUSTM

Linear Cutter





- Uniform compression & consistent staple formation
- Low profile anvil tip to minimise tissue trauma



Cut-line & Staple Line markers for precise cutting & stapling



New knife with every reload guarantees smooth & precise transection



Every reload has a safety lock to prevent firing over empty reload



- Push-Button for quick release
- Ambidextrous-Dual side firing for ease of use




Rear hinge alignment offers ease of handling

Technical Specifications

Stapler	Description	Reloads	Open Length	Close Length	Color	Number of Firings	Number of Staples	Rows of Staple	Staple Line Length	Cut Length
MLC60	Mirus Disp. Linear Cutter 60	MLCR60-3.8 MLCR60-4.8	3.8 mm 4.8 mm	1.5 mm 2.0 mm	Blue Green	8	64	4	64 mm	60 mm
MLC80	Mirus Disp. Linear Cutter 80	MLCR80-3.8 MLCR80-4.8	3.8 mm 4.8 mm	1.5 mm 2.0 mm	Blue Green	8	84	4	84 mm	80 mm
MLC100	Mirus Disp. Linear Cutter 100	MLCR100-3.8 MLCR100-4.8	3.8 mm 4.8 mm	1.5 mm 2.0 mm	Blue Green	8	104	4	104 mm	100 mm

Meril's Global Presence



Meril

Meril Endo-Surgery Pvt. Ltd.
E1-E3, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg, Chala,
Vapi - 396 191, Gujarat, India