



Monik™


Disposable Endoscopic Trocar

Monik™ Disposable Endoscopic Trocar has applications in abdominal, thoracic and gynecological minimally invasive surgical procedures to establish a path of entry for endoscopic instruments. The trocar may be used with or without visualization for primary and secondary insertions.



**Optical Entry :**

Enables visualization of tissue layers during insertion

**Bird Wing Tip :**

Separates the tissue rather than cutting it

**Inner & Outer Seal :**

Maintains insufflation & accepts a wide range of instruments sizes

**Stability Thread :**

Abdominal wall retention for minimal trocar slip-outs

**Obturator :**

Robust design for insertion in abdominal wall with less force

Technical Specifications

Stapler	Description	Sleeve Length
MNK-5	Bladeless Trocar, 5 mm diameter	95 mm
MNK-10	Bladeless Trocar, 10 mm diameter	95 mm
MNK-12	Bladeless Trocar, 12 mm diameter	95 mm
MNK-15	Bladeless Trocar, 15 mm diameter	110 mm