




Appropriate tissue capture and penetration for secure fixation at multiple angles

PROFOUND™ A

Poly (lactide-co-glycolide)




Tube Diameter $5.50 \pm 0.2\text{mm}$



1. Trust, Delivered

Trusted from its design, to its delivery, to its easy to see color.



Screw design spikes with sharp exposed tip ensures minimal tissue damage.

Wide proximal wings holds larger pore size mesh as well as dual meshes securely, ensuring strong repair.




2. Optimal Visibility, Delivered
Violet spike



- Designed specifically to optimize visibility during surgery

3. More Value, Delivered
More Spikes, More Choices



Product Specifications

Profound™ A - Absorbable Mesh Fixation Device

Sr. No.	Product Code	Brand Name	Description
1	MFD 15A	PROFOUND™ A	PLGA - ABSORBABLE (15)
2	MFD 30A	PROFOUND™ A	PLGA - ABSORBABLE (30)